

Product datasheet for SC334565

ANKRD10 (NM_001286721) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ANKRD10 (NM_001286721) Human Untagged Clone
Tag:	Tag Free
Symbol:	ANKRD10
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>NCBI ORF sequence for NM_001286721, the custom clone sequence may differ by one or more nucleotides

ATGTCGGCGGCGGGAGCGGGCGCGGGCGTAGAGGCGGGCTTCTCCAGCGAGGAGCTGCTCTCGCTCCGTT
 TCCCGCTGCACCGCGCCTGCCGCGACGGGGACCTGGCCACGCTCTGCTCGCTGCTGCAGCAGACACCCCA
 CGCCACCTGGCCTCTGAGGACTCCTTCTATGGCTGGACGCCCGTGCACTGGGCCGCGCATTTTCGGCAAG
 TTGGAGTGCTTAGTGCACTTGGTGAGAGCGGGAGCCACACTCAACGCTCTCCACCACCGGTACGCGCAGA
 CGCCAGCCACATTGCAGCCTTTGGGGGACATCCTCAGTGCCTGGTCTGGCTGATTCAAGCAGGAGCCAA
 CATTAAACAAACCGGATTGTGAGGGTGAACTCCATTACAAGGCAGCTCGCTCTGGGAGCCTAGAATGC
 ATCAGTGCCCTTGTGGCAATGGGGCTCAGTCGAGTTTATTTACAGAGCAGTGTACATTCTTGCTTTC
 CAGGGGAACCTTCAACATGGAGTTACTTTTGATCCCTCAGTTTTAATTCAGTGTCTAAAGCCTGAGAAATG
 CCAGTGGCCTGACAGCAGCAGACATTGCACAAACCCAGGGTTTCCAAGAGTGTGCCAGTTTCTCTTGAA
 CCTCCAGAATTGTCATCTGAACCATTCTATAA

Restriction Sites:	SgfI-MluI
ACCN:	NM_001286721
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).


[View online »](#)

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: NM_001286721.1, NP_001273650.1

RefSeq Size: 1310 bp

RefSeq ORF: 663 bp

Locus ID: 55608

UniProt ID: Q9NXR5

Cytogenetics: 13q34