

Product datasheet for **SC334577**

TTC39A (NM_001297666) Human Untagged Clone

Product data:

| | |
|---------------------------|--|
| Product Type: | Expression Plasmids |
| Product Name: | TTC39A (NM_001297666) Human Untagged Clone |
| Tag: | Tag Free |
| Symbol: | TTC39A |
| Synonyms: | C1orf34; DEME-6 |
| Mammalian Cell Selection: | Neomycin |
| Vector: | pCMV6-Entry (PS100001) |
| E. coli Selection: | Kanamycin (25 ug/mL) |
| Fully Sequenced ORF: | >NCBI ORF sequence for NM_001297666, the custom clone sequence may differ by one or more nucleotides |

```
ATGAAGGCCGCTACCTCAGCATGTTTGGGAAGGAGACCACAAGCCGTTCTGGGGACGACGAAGTGAAT
TATTTTCGAGCTGTGCCAGGCCTGAAGCTCAAGATTGCTGGGAAATCTCTACCCACAGAGAAGTTTGCCT
CCGGAAGTCCCGGCGCTACTTCTCCTCAAACCTATCTCGCTGCCAGTGCCTGCTCTGGAATGATGTAC
ATCTGGAACGGCTACGCCGTGATTGGGAAGCAGCCGAACTCACGGATGGGATACTTGAGATTATCACTA
AGGCTGAAGAGATGCTGGAGAAAGGCCAGAGAACGAGTACTCAGTGGATGACGAGTCTTGGTGAATTT
GTTGAAAGGCCTGTCTGAAATACCTGGGCCGTGTCCAGGAGCCGAGGAGAATTTTAGGAGCATCTCT
GCCAATGAAAAGAAGATTAATATGACCACTACTTGATCCCAAACGCCCTGCTGGAGCTGGCCCTGCTGC
TTATGGAGCAAGACAGAAACGAAGAGGCCATCAAACCTTTTGGAACTCTGCCAAGCAAACTACAAGAATTA
CTCCATGGAGTCAAGGACACACTTTCGAATCCAGGCAGCCACACTCCAAGCCAAGTCTTCCCTAGAGAAC
AGCAGCAGATCCATGGTCTCATCAGTGTCTTAG
```

| | |
|--------------------|--|
| Restriction Sites: | Sgfl-MluI |
| ACCN: | NM_001297666 |
| OTI Disclaimer: | Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP). |
| Components: | The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water). |



[View online »](#)

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001297666.1](#), [NP_001284595.1](#)

RefSeq Size: 1769 bp

RefSeq ORF: 666 bp

Locus ID: 22996

UniProt ID: [Q5SRH9](#)

Cytogenetics: 1p32.3