

Product datasheet for **SC334608**

SURF4 (NM_001280790) Human Untagged Clone

Product data:

| | |
|----------------------------------|--|
| Product Type: | Expression Plasmids |
| Product Name: | SURF4 (NM_001280790) Human Untagged Clone |
| Tag: | Tag Free |
| Symbol: | SURF4 |
| Synonyms: | ERV29 |
| Mammalian Cell Selection: | Neomycin |
| Vector: | pCMV6-Entry (PS100001) |
| E. coli Selection: | Kanamycin (25 ug/mL) |
| Fully Sequenced ORF: | >NCBI ORF sequence for NM_001280790, the custom clone sequence may differ by one or more nucleotides |

```
ATGTGGTTCCAGTGGAGCGAGCAGCGCGACTACATCGACACCACCTGGAAGTGGGCTACCTGCTGGCCT
CGTCTTCCTGCTTCTCAACTTGTGGGACAGCTGACTGGCTGCGTCTGGTGTGAGCAGGAACCTCGT
GCAGTACGCCCTGCTTCGGGCTCTTTGGAATCATAGCTCTGCAGACGATTGCCTACAGCATTTTATGGGAC
TTGAAGTTTTTATGAGGAACCTGGCCCTGGGAGGAGGCCTGTTGCTGCTCTAGCAGAATCCCCTTCTG
AAGGGAAGAGCATGTTTGCGGCGTCCCCACCATGCGTGAGAGCTCCCCAAACAGTACATGCAGCTCGG
AGGCAGGGTCTTGCTGGTCTGATGTTTCATGACCTCCTTCACTTTGACGCCAGCTTCTTTCTATTGTC
CAGAACATCGTGGGCACAGCTCTGATGATTTTAGTGGCCATTGGTTTTAAACCAAGCTGGCTGCTTTGA
CTCTTGTGTGTGGCTCTTTGCCATCAACGTATATTTCAACGCCTTCTGGACCATTCCAGTCTACAAGCC
CATGCATGACTTCTGAAATACGACTTCTCCAGACCATGTCGGTGATTGGGGGCTTCTCTGGTGGTG
GCCCTGGGCCCTGGGGGTGTCTCCATGGATGAGAAGAAGAAGGAGTGTTAA
```

| | |
|---------------------------|--|
| Restriction Sites: | Sgfl-MluI |
| ACCN: | NM_001280790 |
| OTI Disclaimer: | Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP). |
| Components: | The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water). |



| | |
|-------------------------------|---|
| Reconstitution Method: | <ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C. |
| RefSeq: | <u>NM_001280790.1, NP_001267719.1</u> |
| RefSeq Size: | 3218 bp |
| RefSeq ORF: | 681 bp |
| Locus ID: | 6836 |
| Cytogenetics: | 9q34.2 |
| Protein Families: | Transmembrane |
| Gene Summary: | <p>This gene is located in the surfeit gene cluster, which is comprised of very tightly linked housekeeping genes that do not share sequence similarity. The encoded protein is a conserved integral membrane protein that interacts with endoplasmic reticulum-Golgi intermediate compartment proteins. Disruption of this gene results in reduced numbers of endoplasmic reticulum-Golgi intermediate compartment clusters and redistribution of coat protein I to the cytosol. Alternate splicing results in multiple transcript variants. [provided by RefSeq, Jul 2013]</p> <p>Transcript Variant: This variant (4) contains an alternate exon in its 5' UTR and initiates translation at a downstream start codon, compared to variant 1. The encoded isoform (4) has a shorter N-terminus, compared to isoform 1. Both variants 4 and 5 encode the same isoform.</p> |