

Product datasheet for SC334653

DCAKD (NM_001288655) Human Untagged Clone

Product data:

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

| Product Type: | Expression Plasmids |
|---------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name: | DCAKD (NM_001288655) Human Untagged Clone |
| Tag: | Tag Free |
| Symbol: | DCAKD |
| Vector: | pCMV6 series |
| Fully Sequenced ORF: | >NCBI ORF sequence for NM_001288655, the custom clone sequence may differ by one or more nucleotides |
| | ATGTTTCTGGTGGGCCTGACAGGGGGCATTGCCTCAGGCAAGAGCTCAGTGATCCAGGTGTTCCAGCAGC TGGGCTGTGCGGTGATTGACGTGGACGTGATGGCCCGGCACGTCGTGCAGCCAGGATACCCTGCCACCG GCGCATCGTAGAGGTCTTCGGCACTGAGGTCTTGCTGGAGAACGGCGACATAAATCGCAAGGTCCTGGGG GACCTGATCTTTAACCAGCCTGACCGGCGGCAGCTGCTCAACGCCATCACCCACC |
| Restriction Sites: | Sgfl-Mlul |
| ACCN: | NM_001288655 |
| OTI Disclaimer: | Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP). |
| Components: | The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water). |



This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2022 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US

DCAKD (NM_001288655) Human Untagged Clone – SC334653

| Reconstitution Method: | Centrifuge at 5,000xg for 5min. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. Close the tube and incubate for 10 minutes at room temperature. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C. |
|------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| RefSeq: | <u>NM 001288655.1, NP 001275584.1</u> |
| RefSeq Size: | 2072 bp |
| RefSeq ORF: | 696 bp |
| Locus ID: | 79877 |
| UniProt ID: | <u>Q8WVC6</u> |
| Cytogenetics: | 17q21.31 |
| Protein Families: | Transmembrane |

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2022 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US