

## Product datasheet for **SC334670**

### PAQR6 (NM\_001272108) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	PAQR6 (NM_001272108) Human Untagged Clone
Tag:	Tag Free
Symbol:	PAQR6
Synonyms:	PRdelta
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>NCBI ORF sequence for NM_001272108, the custom clone sequence may differ by one or more nucleotides

```
ATGCTCAGTCTCAAGCTGCCCAACTTCTCAAGTCCACCAGGTCCCCGGGTGTTCTGGGAAGATGGCA
TCATGTCTGGCTACCGCCGCCACCAGCTCGGCTTTGGACTGTGTCTCAGCTCCTCCAGATGACCAA
CGAGACGGTCAACATCTGGACTCACTTCTGCCACCTGTTTCTGGAGCTGAAAGCCCTGGGCTCAGT
AAGGTCTCCGCACAGGAGCCTTCGCCTATCCATTCTGTTCGACAACCTCCCCTCTTTTATCGGCTCG
GGCTGTGCTGGGGCAGGGGCCACGGCTGTGGGCAGGAGGCCCTGAGCACCAGCCATGGCTACCATCTTT
CTGCGCGCTGCTCACTGGCTTCTTTCGCCTCCCACCTGCCTGAAAGGCTGGCACCAGGACGCTTTGAT
TACATCGGCCACAGCCACCAGTTATTCCACATCTGTGCAGTGCTGGGCACCCACTCCAGCTGGAGGCAG
TGCTGGCTGATATGGGATCACGCAGAGCCTGGCTGGCCACACAGGAACCTGCCCTGGGCCTGGCAGGCAC
AGTGGCCACACTGGTCTTTGGCTGCAGCTGGGAACCTACTCATTATTGCTGCTTTACAGCCACCCTGCTT
CGGGCCCCAGTACATGCCCTCTGCTGCAGGGTGGCCACTGGAGGGGGGTACCCAGGCCAAACAACAGT
GA
```

Restriction Sites:	Sgfl-MluI
ACCN:	NM_001272108
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



[View online »](#)

<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u><a href="#">NM_001272108.1</a></u> , <u><a href="#">NP_001259037.1</a></u>
<b>RefSeq Size:</b>	1696 bp
<b>RefSeq ORF:</b>	702 bp
<b>Locus ID:</b>	79957
<b>UniProt ID:</b>	<u><a href="#">Q6TCH4</a></u>
<b>Cytogenetics:</b>	1q22
<b>Protein Families:</b>	Druggable Genome, Transmembrane
<b>Gene Summary:</b>	<p>Plasma membrane progesterone (P4) receptor coupled to G proteins (PubMed:23763432, PubMed:23161870). Seems to act through a G(s) mediated pathway (PubMed:23161870). Involved in neurosteroid inhibition of apoptosis (PubMed:23161870). May be involved in regulating rapid P4 signaling in the nervous system (PubMed:23763432). Also binds dehydroepiandrosterone (DHEA), pregnanolone, pregnenolone and allopregnanolone (PubMed:23763432, PubMed:23161870).[UniProtKB/Swiss-Prot Function]</p> <p>Transcript Variant: This variant (7) has multiple differences and initiates translation at an alternate start codon, compared to variant 1. The resulting isoform (6) has distinct N- and C-termini and is shorter, compared to isoform 1.</p>