

Product datasheet for **SC334681**

TMEM249 (NM_001280561) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	TMEM249 (NM_001280561) Human Untagged Clone
Tag:	Tag Free
Symbol:	TMEM249
Synonyms:	C8ORFK29
Vector:	<u>pCMV6 series</u>
Fully Sequenced ORF:	>NCBI ORF sequence for NM_001280561, the custom clone sequence may differ by one or more nucleotides

ATGCCTAAGGGCCGGCCGGCAGCCTGCCACGACCTCCATCGGCTGGCGCTTCCAGCTCTGGTTCTAG
 GACTCACCTGCCCGAGCGCCACCTGGCCAGACGCCTCAAGAACAATAGCTTCTACCCATTCGTGCAGCA
 GGAGCCAAACGTCTTCGTGCTCGAATACTACCTGGACACACTGTGAAGGGGATGCTGCTTTCATCATC
 AGCGTGGTCCTGGTCAGCTTCAGCTCCCTGAGAGAGGTGCAGAAGCAAGAGACCTGGGTCTTCCTCGTGT
 ACGGCGTGGGCGTGGGCTGTGGCTCGTGATCTCGTCGCTGCCGCGGCGCCGCTGGTGTGAACACAC
 GCGCGGCGTGTACCACTTCTCCATCCAGGCGCGACCGTGTGCCAGGGCCCTTGATCTGGTCTACGTG
 CGCCTGGCGCTCAGCTCCGATGCCACGGAAGGTGCTTCTCCACCTGGTCTCGGCGGCCACAGTTGG
 AGCCCCCTGGTGTGGTGCAGCTGTGGAGCACTACGAGCAAAATGGAATACTTGGGCCGTTACATCGCCCC
 GAAACTCAACATCAACTACTTCGACTACCTGGCCACCTCCTACCGGCACGTGGTTCGCCACTGGCCACCC
 CCTGGCGCTGGCACTGTGATGGCAAGAGCCCCATGGGTCACAAGCCCAGCTCCTCTCAGTCCAGCCTGG
 AGGTGTGA

Restriction Sites:	SgfI-MluI
ACCN:	NM_001280561
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).


[View online »](#)

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: NM_001280561.1, NP_001267490.1

RefSeq Size: 1598 bp

RefSeq ORF: 708 bp

Locus ID: 340393

UniProt ID: Q2WGJ8

Cytogenetics: 8q24.3