

Product datasheet for **SC334754**

KCNIP1 (NM_001278339) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	KCNIP1 (NM_001278339) Human Untagged Clone
Tag:	Tag Free
Symbol:	KCNIP1
Synonyms:	KCHIP1; VABP
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>NCBI ORF sequence for NM_001278339, the custom clone sequence may differ by one or more nucleotides

```
ATGGGGCCGTCATGGGCACCTTCTCATCTCTGCAAACCAACAAAGGCGACCCTCGAAAGATAAGATTG
AAGATGAGCTGGAGATGACCATGGTTTGCATCGGCCCGAGGGACTGGAGCAGCTCGAGGCCAGACCAA
CTTCACCAAGAGGGAGCTGCAGGTCCTTTATCGAGGCTTCAAAAATGAGTGCCCCAGTGGTGTGGTCAAC
GAAGACACATTCAAGCAGATCTATGCTCAGTTTTCCCTCATGGAGCCCTGCCTTGTGGAAAGTTCTC
CCTGTGTGGAATTCCTGCCCCATCACCTGCCCTCCTTTTCTGCCTTGTAGATGCCAGCACGTATGCCCA
TTACCTCTTCAATGCCTTCGACACCACTCAGACAGGCTCCGTGAAGTTCGAGGACTTTGTAACCGCTCTG
TCGATTTTATTGAGAGGAAGTCCACGAGAACTAAGGTGGACATTTAATTTGTATGACATCAACAAGG
ACGGATACATAAACAAGAGGAGATGATGGACATTGTCAAAGCCATCTATGACATGATGGGAAATACAC
ATATCCTGTGCTCAAAGAGGACTCCAAGGCAGCATGTGGACGTCTTCTCCAGAAAATGGACAAAAAT
AAAGATGGCATCGTAACTTTAGATGAATTTCTGAATCATGTCAGGAGGACGACAACATCATGAGGTCTC
TCCAGCTGTTTCAAAATGTCATGTAA
```

Restriction Sites:	Sgfl-MluI
ACCN:	NM_001278339
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



[View online »](#)

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001278339.1</u> , <u>NP_001265268.1</u>
RefSeq Size:	2103 bp
RefSeq ORF:	726 bp
Locus ID:	30820
UniProt ID:	<u>Q9NZI2</u>
Cytogenetics:	5q35.1
Protein Families:	Druggable Genome, Ion Channels: Other
Gene Summary:	<p>This gene encodes a member of the family of cytosolic voltage-gated potassium (Kv) channel-interacting proteins (KCNIPs), which belong to the neuronal calcium sensor (NCS) family of the calcium binding EF-hand proteins. They associate with Kv4 alpha subunits to form native Kv4 channel complexes. The encoded protein may regulate rapidly inactivating (A-type) currents, and hence neuronal membrane excitability, in response to changes in the concentration of intracellular calcium. Alternative splicing results in multiple transcript variants encoding different isoforms. [provided by RefSeq, May 2013]</p> <p>Transcript Variant: This variant (5) encodes the longest isoform (5).</p>