

Product datasheet for **SC334825**

SERHL2 (NM_001284334) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	SERHL2 (NM_001284334) Human Untagged Clone
Tag:	Tag Free
Symbol:	SERHL2
Synonyms:	dJ222E13.1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>NCBI ORF sequence for NM_001284334, the custom clone sequence may differ by one or more nucleotides

ATGAGTGAGAACGCCGACCAGACTTTTATTACGTTGCCATGGATTTTCGGAGGTCATGGGCTCTCGTCCC
 ATTACAGCCCAGGTGTCCCATATTACCTCCAGACTTTTGTGAGTGAGATCCGAAGAGTTGTGGCAGGTGG
 CGTCGTGGGCGGAATGTTTTCTGTACCTTCCCCGAGATGGTGGATAAACTTATCTTGCTGGACACGCCG
 CTCTTTCTCCTGGAATCAGATGAAATGGAGAACTTGCTGACCTACAAGCGGAGAGCCATAGAGCACGTGC
 TGCAGGTAGAGGCCTCCCAGGAGCCCTCGCACGTGTTGAGCCTGAAGCAGCTGCTGCAGAGGTTACTGAA
 GAGCAATAGCCACTTGAGTGAGGAGTGCGGGGAGCTTCTCCTGCAAAGAGGAACCAAGGTGGCCACA
 GGTCTGGTTCTGAACAGAGACCAGAGGCTCGCCTGGGCAGAGAACAGCATTGACTTCATCAGCAGGGAGC
 TGTGTGCGCATTCCATCAGGAAGCTGCAGGCCCATGTCTGTTGATCAAAGCAGTCCACGGATATTTTGA
 TTCAAGACAGAATTACTCTGAGAAGGAGTCCCTGTCGTTTCATGATAGACACGATGAAATCCACCCTCAAA
 GAGCAGTTCAGTTTGTGGAAGTCCCAGGCAATCACTGTGTCCACATGAGCGAACCCAGCACGTGGCCA
 GTATCATCAGCTCCTTCTTACAGTGCACACACATGCTCCAGCCAGCTGTAG

Restriction Sites:	Sgfl-MluI
ACCN:	NM_001284334
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).


[View online »](#)

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001284334.1, NP_001271263.1</u>
RefSeq Size:	1255 bp
RefSeq ORF:	753 bp
Locus ID:	253190
UniProt ID:	<u>Q9H4I8</u>
Cytogenetics:	22q13.2
Gene Summary:	Probable serine hydrolase. May be related to cell muscle hypertrophy.[UniProtKB/Swiss-Prot Function] Transcript Variant: This variant (2) lacks two alternate in-frame exons in the 5' coding region, compared to variant 1, resulting in an isoform (b) that is shorter than isoform a.