

## Product datasheet for **SC334902**

### **IKZF3 (NM\_001284514) Human Untagged Clone**

#### Product data:

Product Type:	Expression Plasmids
Product Name:	IKZF3 (NM_001284514) Human Untagged Clone
Tag:	Tag Free
Symbol:	IKZF3
Synonyms:	AIO; AIOLOS; ZNFN1A3
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>NCBI ORF sequence for NM_001284514, the custom clone sequence may differ by one or more nucleotides

```
ATGGGAAGTCAAAGAGCTCTCGTACTGGACAGATTAGCAAGCAATGTGGCAAAACGAAAAAGCTCAATGC
CTCAGAAATTCATTGGTGAAGCGCCACTGCTTTGATGTCAACTATAATTCAAGTTACATGTATGAGAA
AGAGAGTGAGCTCATAACAGACCCGCATGATGGACCAAGCCATCAATAACGCCATCAGCTATCTTGGCGCC
GAAGCCCTGCGCCCTTGGTCCAGACACCCGCTGCTCCCACCTCGGAGATGGTCCAGTTATCAGCAGCA
TGTATCCCATAGCCCTCACCCGGCTGAGATGTCAAACGGTGCCCTCAAGAGCTGGAAAAGAAAAGCAT
CCACCTTCCAGAGAAGAGCGTGCCTTCTGAGAGAGGCCTCTCTCCAACAATAGTGGCCACGACTCCAGC
GACTGACAGCAACCATGAAGAAGCCAGAATCACATCTATCAGCAAAATCACATGGTCTGTCTCGGG
CCCGCAATGGGATGCCACTTCTGAAGGAGTTCCCGCTTACGAACTCCTCAAGCCCCGCCCATCTG
CCCAAGAGACTCCGTCAAAGTGATCAACAAGGAAGGGGAGGTGATGGATGTGTATCGGTGTGACCACTGC
CGCGTCTCTTCTGGACTATGTGATGTTACGATTACATGGGCTGCCACGGCTTCCGTGACCCCTTTCG
AGTGTAACATGTGTGGATATCGAAGCCATGATCGGTATGAGTTCTCGTCTCACATAGCCAGAGGAGAACA
CAGAGCCCTGCTGAAGTGA
```

Restriction Sites:	Sgfl-MluI
ACCN:	NM_001284514
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



[View online »](#)

<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u>NM_001284514.1, NP_001271443.1</u>
<b>RefSeq Size:</b>	9812 bp
<b>RefSeq ORF:</b>	789 bp
<b>Locus ID:</b>	22806
<b>UniProt ID:</b>	<u>Q9UKT9</u>
<b>Cytogenetics:</b>	17q12-q21.1
<b>Gene Summary:</b>	<p>This gene encodes a member of the Ikaros family of zinc-finger proteins. Three members of this protein family (Ikaros, Aiolos and Helios) are hematopoietic-specific transcription factors involved in the regulation of lymphocyte development. This gene product is a transcription factor that is important in the regulation of B lymphocyte proliferation and differentiation. Both Ikaros and Aiolos can participate in chromatin remodeling. Regulation of gene expression in B lymphocytes by Aiolos is complex as it appears to require the sequential formation of Ikaros homodimers, Ikaros/Aiolos heterodimers, and Aiolos homodimers. Several alternative transcripts encoding different isoforms have been described, as well as some non-protein coding variants. [provided by RefSeq, Apr 2012]</p> <p>Transcript Variant: This variant (14) uses an alternate 5' exon structure and thus differs in the 5' UTR and 5' coding region, compared to variant 1. These differences cause translation initiation at a downstream start codon and result in an isoform (14) with a shorter N-terminus, compared to isoform 1. Variants 14, 15, and 16 encode the same isoform (14).</p> <p>Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.</p>