

Product datasheet for **SC334947**

PGAP3 (NM_001291726) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	PGAP3 (NM_001291726) Human Untagged Clone
Tag:	Tag Free
Symbol:	PGAP3
Synonyms:	AGLA546; CAB2; hCOS16; PERLD1; PP1498
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>NCBI ORF sequence for NM_001291726, the custom clone sequence may differ by one or more nucleotides

```
ATGGCCGGCCTGGCGGCGCGGTTGGTCCTGCTAGCTGGGGCAGCGGCGCTGGCGAGCGGCTCCCAGGGCG  
ACCGTGAGCCGGTGTACCGGACTGCGTACTGCAAGAGCAGAAGTCTGCTGGGGCGCTCTGAA  
TCACTTCCGCTCCCGCCAGCCAATCTACATGAGTCTAGCAGGCTGGACCTGTCGGGACGACTGTAAGTAT  
GAGTGTATGTGGGTACCGTTGGGCTCTACCTCCAGGAAGGTCAAAAAGTGCCTCAGTTCATGGCAAGG  
TGTCCTCAATGCATGGTCTGGTCCACAGTTTTCCACACCAGGGACTGACCTCACAGAGAAAATGGA  
CTACTTCTGTGCCTCCACTGTCATCTACTCAATCTACCTGTGCTGCGTCAGGACCGTGGGGCTGCAG  
CACCCAGCTGTGGTCAGTGCCTTCGGGCTCTCCTGCTGCTCATGCTGACCGTGCAGTCTCCTACCTGA  
GCCTCATCCGCTTCGACTATGGCTACAACCTGGTGCCAAACGTGGCTATTGGCCTGGTCAACGTGGTGTG  
GTGGCTGGCCTGGTGCCTGTGGAACCAGCGGCGGCTGCCTCACGTGCGCAAGTGCCTGGTGGTGGTCTTG  
CTGCTGCAGGGGCTGTCCCTGCTCGAGCTGCTTACTCCACCCTCTTCTGGGTCTGGATGCCCATG  
CCATCTGGCAGATCAGCACCATCCCTGTCCACGTCCTCTTTTCAGCTTCTGGAAGATGACAGCCTGTA  
CCTGCTGAAGGAATCAGAGGACAAGTTCAAGCTGGACTGA
```

Restriction Sites:	Sgfl-MluI
ACCN:	NM_001291726
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



[View online »](#)

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001291726.1, NP_001278655.1</u>
RefSeq Size:	2568 bp
RefSeq ORF:	810 bp
Locus ID:	93210
UniProt ID:	<u>Q96FM1</u>
Cytogenetics:	17q12
Protein Families:	Transmembrane
Gene Summary:	<p>This gene encodes a glycosylphosphatidylinositol (GPI)-specific phospholipase that primarily localizes to the Golgi apparatus. This ubiquitously expressed gene is predicted to encode a seven-transmembrane protein that removes unsaturated fatty acids from the sn-2 position of GPI. The remodeling of the constituent fatty acids on GPI is thought to be important for the proper association between GPI-anchored proteins and lipid rafts. The tethering of proteins to plasma membranes via posttranslational GPI-anchoring is thought to play a role in protein sorting and trafficking. Mutations in this gene cause an autosomal recessive form of neurologic hyperphosphatasia with cognitive disability (HPMRS4). Alternative splicing results in multiple transcript variants encoding distinct isoforms. [provided by RefSeq, Jul 2017]</p> <p>Transcript Variant: This variant (2) lacks an in-frame exon in the 5' coding region, compared to variant 1. It encodes a shorter protein (isoform 2), compared to isoform 1.</p>