

## Product datasheet for **SC334957**

### PRUNE (PRUNE1) (NM\_001303229) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	PRUNE (PRUNE1) (NM_001303229) Human Untagged Clone
Tag:	Tag Free
Symbol:	PRUNE1
Synonyms:	DRES-17; DRES17; H-PRUNE; HTCD37; NMIHBA; PRUNE
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>NCBI ORF sequence for NM_001303229, the custom clone sequence may differ by one or more nucleotides

```
ATGGACCTTAAAATTGGAAGGCAACCCCAAAGGACAGCAAATATGTGGAGAACTAGAGGCCCTTTTCC
CAGACCTACCCAAGAGAAATGATATATTTGATTCCTACAAAAGGCAAAGTTTGATGATCAGGACTGAC
CACTGAGCAGATGCTGAGAAAAGACCAGAAGACTATCTATAGACAAGGCGTCAAGGTGGCCATTAGTGCA
ATATATATGGATTTGGAGGCCTTTCTGCAGAGGTCTAACCTCCTTGCAGATCTCCATGCTTTCTGCCAGG
CTCACAGCTATGATGTCCTGGTTGCCATGACTATCTTTTTCAACACTCACAATGAGCCAGTGCGGCAGTT
GGCTATTTCTGTCCCCATGTGGCACTCCAACAACGATCTGTGAAGTCTGGAACGCTCCCACTCTCCA
CCCCTGAAGCTGACCCCTGCCTCAAGTACCCACCCTAACCTCCATGCCTATCTTCAAGGCAACCCAGG
TCTCTCGAAAGAACTTCTGCCCTGCTCCAGGAAGCCCTGTGAGCATATTTGACTCCATGAAGATCCC
TTCAGGACAGCCTGAGACAGCAGATGTGTCAGGGAGCAAGTGGACAAGGAATTGGACAGGGCAAGTAAC
TCCCTGATTTCTGGCCTGAGTCAAGATGAGGAGGACCCTCCGCTGCCCCGACGCCCATGAACAGCTTGG
TGGATGAGTGCCTCTAGATCAGGGGCTGCCTAAACTCTCTGCTGAGGCCGTCTTCGAGAAGTGCAGTCA
GATCTCACTGTACAGTCTACCACAGCCTCCCTGTCCAAGAAGTGA
```

Restriction Sites:	Sgfl-MluI
ACCN:	NM_001303229
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



[View online »](#)

<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u>NM_001303229.1, NP_001290158.1</u>
<b>RefSeq Size:</b>	2898 bp
<b>RefSeq ORF:</b>	816 bp
<b>Locus ID:</b>	58497
<b>UniProt ID:</b>	<u>Q86TP1</u>
<b>Cytogenetics:</b>	1q21.3
<b>Protein Pathways:</b>	Purine metabolism
<b>Gene Summary:</b>	<p>This gene encodes a member of the DHH protein superfamily of phosphoesterases. This protein has been found to function as both a nucleotide phosphodiesterase and an exopolyphosphatase. This protein is believed to stimulate cancer progression and metastases through the induction of cell motility. A pseudogene has been identified on chromosome 13. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Dec 2014]</p> <p>Transcript Variant: This variant (2) lacks one alternate internal exon resulting in a distinct 5' UTR and the use of an in-frame downstream start codon compared to variant 1. The resulting protein (isoform 2) has a shorter N-terminus than isoform 1.</p>