

Product datasheet for **SC334969**

PITPNB (NM_001284277) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	PITPNB (NM_001284277) Human Untagged Clone
Tag:	Tag Free
Symbol:	PITPNB
Synonyms:	PI-TP-beta; PtdInsTP; VIB1B
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>NCBI ORF sequence for NM_001284277, the custom clone sequence may differ by one or more nucleotides

```
ATGGTGCTGATCAAGGAATCCGTGTGGTTTTGCCATGTTCTGTTCCAGGAGTATCAGGTTGGGCAGCTTT
ACTCTGTTGCAGAAGCTAGTAAGAATGAGACTGGTGGTGGAGAAGGAATTGAAGTCTTAAAGAATGAACC
TTATGAGAAGGATGGAGAAAAGGACAGTATACGCACAAAATTTATCACCTAAAGAGCAAAGTGCCTGCA
TTCGTGAGGATGATTGCTCCCGAGGGCTCCTTGGTGTTCATGAGAAAGCCTGGAATGCGTACCCCTACT
GTAGAACAATTGTAAACGAATGAATATATGAAAGATGATTTCTTCATTAATAATCGAAACATGGCACAACC
AGACTTGGGAACATTAGAAAATGTACATGGTTTAGATCCAAACACATGGAAAACCTGTTGAAATTGCCAT
ATAGATATTGCAGATAGAAGTCAAGTTGAACCAGCAGACTACAAAGCTGATGAAGACCCAGCATTATTCC
AGTCAGTCAAGACCAAGAGAGGCCCTTTGGGACCAACTGGAAGAAGGAGCTGGCAAACAGCCCTGACTG
TCCCCAGATGTGTGCCTATAAGCTGGTGACCATCAAATCAAGTGGTGGGACTGCAAAGCAAAGTAGAA
AACTTCATTCAAAGCAAGAAAAACGGATATTTACAAACTTCCATCGCCAGCTTTTTTGTGGATTGACA
AGTGGATCGATCTCAGTGAAGACATTAGGAGAATGGAAGACGAGACTCAGAAAGAACTAGAAACATT
ACGTAATCAAGGTCAAGTGAGGGGAACAAGCGCAGCCAGTATGAGTGA
```

Restriction Sites:	Sgfl-MluI
ACCN:	NM_001284277
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



[View online »](#)

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001284277.1, NP_001271206.1</u>
RefSeq Size:	2895 bp
RefSeq ORF:	819 bp
Locus ID:	23760
UniProt ID:	<u>P48739</u>
Cytogenetics:	22q12.1
Gene Summary:	<p>This gene encodes a cytoplasmic protein that catalyzes the transfer of phosphatidylinositol and phosphatidylcholine between membranes. This transfer activity is required for COPI complex-mediated retrograde transport from the Golgi apparatus to the endoplasmic reticulum. Alternative splicing of this gene results in multiple transcript variants. [provided by RefSeq, Sep 2013]</p> <p>Transcript Variant: This variant (sp2) lacks an exon that results in an alternate 3' coding region, compared to variant sp1. The encoded isoform (2) has a distinct C-terminus and is 1 aa longer than isoform 1.</p>