

Product datasheet for **SC335023**

PPIL6 (NM_001286360) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	PPIL6 (NM_001286360) Human Untagged Clone
Tag:	Tag Free
Symbol:	PPIL6
Synonyms:	bA425D10.6; dj919F19.1; PPIase; RSPH12
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>NCBI ORF sequence for NM_001286360, the custom clone sequence may differ by one or more nucleotides

```
ATGGCAAGGCCGACGCCGTGCGGGCCCCGCACGCTAGGTGCGGCTCGCCGTCGCTGCCGGAGCGGCCGC
TGCAGGTGAAGGTGGTGGGGCTCTTCAGCTGCCCAACTTTCAGATTGCGAAGAGCGCCGCTGAGGAACT
CAAAAATGAAACCTGGGAATATTCTTCTCTGTGATTTCTTTTGTAAATGGTCAGTTTCTGGGTGATGCA
TTGGATCTGCAGAAATGGGCCACGAGGTGTGGGATATAGTTGACATAAACCCCTCTGCACTTTATGACG
CACTCACTGAGGATTTTCCGCTAAGTTCTTAAGAGACACCAAGCATGATTCGTGTTTTGGACATTTG
TATTGATCTTCTCCAATTGGAAGATTGATTTTTGAGCTATACTGTGATGTGTGCCAAAACATGTAAA
AATTTTCAGGTCTTGTGCACAGGAAAAGCAGGGTTTTCTCAACGTGGCATAAGACTACATTACAAAATT
CCATTTTTCATCGAATAGTACAGAATGGCTGGATACAAGGAGGGGATATAGTGTATGGAAAAGGAGATAA
TGGAGAGTCGATTTATGGTCCAACATTTGAAGATGAAAACTTTTCAAGTTCCATATAATAAAGAGGAGTA
CTTGGAATGGCCAACAAAGGCCGTCACAGCAACGGGTCACAATTCATATCACACTGCAAGCAACTCCTT
ATCTAGATAGAAAATTTGTGGCTTTTGGGCAACTGATTGAAGGAACAGAAGTGCTTAAACAACACTAGAATT
AGTTCCAACACAGAATGAAAGACCAATACATATGTGTAGAATTACTGACAGTGGAGATCCTTATGCTGA
```

Restriction Sites:	Sgfl-MluI
ACCN:	NM_001286360
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



[View online »](#)

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001286360.1, NP_001273289.1</u>
RefSeq Size:	4049 bp
RefSeq ORF:	840 bp
Locus ID:	285755
UniProt ID:	<u>Q8IXY8</u>
Cytogenetics:	6q21
Gene Summary:	<p>Probable inactive PPlase with no peptidyl-prolyl cis-trans isomerase activity.[UniProtKB/Swiss-Prot Function]</p> <p>Transcript Variant: This variant (3) lacks an alternate in-frame exon in the 5' coding region, compared to variant 1, resulting in an isoform (3) that is shorter than isoform 1. Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.</p>