

Product datasheet for **SC335028**

METT5D1 (METTL15) (NM_001297775) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	METT5D1 (METTL15) (NM_001297775) Human Untagged Clone
Tag:	Tag Free
Symbol:	METT5D1
Synonyms:	METT5D1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>NCBI ORF sequence for NM_001297775, the custom clone sequence may differ by one or more nucleotides

```
ATGCTTCGGTATCCATATTTTTGTAGAATGTATAAAGAATGCCTTTCATGTTGGTTGGAATCTGGCATACTAATTTAGGTGTCTGGCCAAACAGAATACATACTACAGCAGAAAAATATAGAGAATATGAAGCCCGGAGCAAAACAGATCAAACCTCAAGCCAGGAGTTACACAGATCTCAAGATAGAGATTTTGAAACTATGGCTAAATTACATATTCCAGTAATGGTGGATGAAGTTGTTTCATTGTTTGTCCACCACAAAAGGACAGATTTTCTAGATATGACATTTGGTTCGGGAGGGCACACAAAAGCCATTCTGCAGAAGGAGTCAGATATTGTTCTCTATGCTTGGACAGAGACCCAACAGCTTATGCATTAGCTGAACATCTTTCAGAGTTGTATCCTAAACAAATCCGAGCTATGCTGGGCCAGTTCAGCCAGGCAGAAGCCTTATTAATGAAAGCTGGAGTGCAGCCAGGAACTTTTGATGGAGTTCCTATGGATCTTGGGTGTTCCCTCCATGCAACTTGATACTCCTGAAAGAGTTTTCCCTTCGAAAGATGGCCCTTTGGACATGAGAATGGATGGTGGCAGATCAACAGGCACTTGCATCTATCCTAAGAACATACGGGGAGGAGAAGCATGCCAAGAAAATCGCTTCAGCAATTGTTTCAGGCACGCAGCATCTACCCCATCACCAGAACCAGCAGCTTGCCAGCATCGTTGCAGGAGCATTTCCCTCCCTGCTATTTATACACGGAAAAGACTTACTACAGCGATCTACCCATATTGCCACCAAGACTTCCAGGCTCTTCGCATATTTGTGAACAATGA
```

Restriction Sites:	Sgfl-MluI
ACCN:	NM_001297775
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



[View online »](#)

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001297775.1, NP_001284704.1</u>
RefSeq Size:	4301 bp
RefSeq ORF:	840 bp
Locus ID:	196074
UniProt ID:	<u>A6NJ78</u>
Cytogenetics:	11p14.1
Gene Summary:	<p>Probable S-adenosyl-L-methionine-dependent methyltransferase.[UniProtKB/Swiss-Prot Function]</p> <p>Transcript Variant: This variant (3) uses an alternate in-frame splice junction at the 5' end of an exon compared to variant 1. The resulting isoform (3) has a shorter and distinct C-terminus compared to isoform 1. Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.</p>