

## Product datasheet for **SC335068**

### FAM82B (RMDN1) (NM\_001286719) Human Untagged Clone

#### Product data:

|                           |                                                                                                      |
|---------------------------|------------------------------------------------------------------------------------------------------|
| Product Type:             | Expression Plasmids                                                                                  |
| Product Name:             | FAM82B (RMDN1) (NM_001286719) Human Untagged Clone                                                   |
| Tag:                      | Tag Free                                                                                             |
| Symbol:                   | FAM82B                                                                                               |
| Synonyms:                 | CGI-90; FAM82B; RMD-1; RMD1                                                                          |
| Mammalian Cell Selection: | Neomycin                                                                                             |
| Vector:                   | pCMV6-Entry (PS100001)                                                                               |
| E. coli Selection:        | Kanamycin (25 ug/mL)                                                                                 |
| Fully Sequenced ORF:      | >NCBI ORF sequence for NM_001286719, the custom clone sequence may differ by one or more nucleotides |

ATGGCGCTGGCTGCTCGACTGTGGCGCCTTCTGCCTTTCGACGTGGAGCCGCCCGGGGTCTCGTCTCC  
 CTGCGGGGACTTCGGGCAGCCGCGGGCATTGCGGCCCTGTGATTCCGCGGCTTCGAGGTAATGGGAAA  
 CCCAGGAACTTTCAAAGAGGCCCTTTACTCTCAGCTTTGTCGATTTGGGTTTTGAACTTACCAGGTT  
 ATCTCTCAGGCTGCTGTGGTTCATGCCACAGCCAAAGTTGAAGAAATACTGAACAAGCAGACTACCTGT  
 ATGAAAGCGGAGAAACAGAAAACTTTATCAGTTGCTAACCCAATACAAGGAAAGTGAAGATGCAGAGTT  
 ACTGTGGCGTTTGGCACGGGCATCACGTGATGTAGCTCAGCTTAGCAGAACCTCAGAAGAGGAGAAAAAG  
 CTATTGGTGTATGAAGCCCTAGAGTATGCAAAAAGAGCACTAGAAAAAATGAATCAAGTTTGCATCTC  
 ATAAGAAAGCAATTGAACTGAACCCTAAAGATGCTACTTCAATTCACCTTATGGGTATTTGGTGCTATAC  
 ATTTGCCGAAATGCCTTGGTATCAAAGAAGATTGCTAAAATGCTGTTTGCAACTCCTCCTAGTTCCACC  
 TATGAGAAGGCCTTAGGCTACTTTCACAGGGCAGAACAAAGTGGATCCAACTTCTACAGCAAAACTTAC  
 TTCTTTTAGGAAAGACATACTTGAACTACACAACAAAAGCTTGCTGCTTTCTGGCTAATGAAAGCCAA  
 GGACTATCCAGCACACAGAGGAGGATAAACAGATACAGACAGAAGCTGCTCAGTTGCTTACAAGTTTC  
 AGTGAGAAGAATTGA

|                    |                                                                                                                                                                                                                                                                                                                                                |
|--------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Restriction Sites: | SgfI-MluI                                                                                                                                                                                                                                                                                                                                      |
| ACCN:              | NM_001286719                                                                                                                                                                                                                                                                                                                                   |
| OTI Disclaimer:    | Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP). |


[View online »](#)

|                               |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
|-------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Components:</b>            | The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).                                                                                                                                                                                                                                                                                                                              |
| <b>Reconstitution Method:</b> | <ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol> |
| <b>RefSeq:</b>                | <u><a href="#">NM_001286719.1</a></u> , <u><a href="#">NP_001273648.1</a></u>                                                                                                                                                                                                                                                                                                                                                                                                                                     |
| <b>RefSeq Size:</b>           | 2991 bp                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| <b>RefSeq ORF:</b>            | 855 bp                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| <b>Locus ID:</b>              | 51115                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
| <b>UniProt ID:</b>            | <u><a href="#">Q96DB5</a></u>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                     |
| <b>Cytogenetics:</b>          | 8q21.3                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                            |