

Product datasheet for **SC335104**

RFPL4AL1 (NM_001277397) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	RFPL4AL1 (NM_001277397) Human Untagged Clone
Tag:	Tag Free
Symbol:	RFPL4AL1
Synonyms:	RNF210B
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>NCBI ORF sequence for NM_001277397, the custom clone sequence may differ by one or more nucleotides

ATGGCTGAGCACTTCAAACAGATCATTAGATGTCCTGTCTGTCTAAAAGATCTTGAAGAAGCCGTGCAAC
 TGAAATGTGGATATGCCTGCTGCCTCCAGTGCCTCAATCACTCCAGAAGGAGCCGATGGGGAAGGTTT
 ACTGTGCCGTTTCTGCTCTGTGGTCTCTCAGAAGGATGACATCAAGCCCAAGTACAAGCTGAGGGCGCTG
 GTTTCATCATCAAGGAAGTGAAGCCCAAGCTGAAATCTGTTCTAACAATGAACCAAGGATGAGGAAGT
 TTCAAGTGGATATGACCTTCGATGTGGACACAGCCAACAATCTCATCATTTCTGAAGACCTGAGGAG
 TTTCCGAAGTGGGGATTTGAGCCAGAATAGGAAGGAGCAAGCTGAGAGGTTGACACTGCCCTGTGCGTC
 CTGGGCACCCCTCGCTTCACTTCCGGCCGCCATTACTGGGAGGTGGACGTGGGCACCAGCCAAGTGTGGG
 ATGTGGGCGTGTGCAAGGAATCTGTCAACCGACAGGGGAAGATTGAGCTTTCTTCAGAACACGGCTTCTT
 GACTGTGGGTTGCAGAGAAGGAAGGTCTTTGCTGCCAGCACTGTGCCTATGACTCCTCTCTGGGTGAGT
 CCCCAGTTGCACAGAGTGGGGATTTTCCTGGATGTAGGTATGAGGTCCATTGCCTTTTACAATGTTAGTG
 ATGGGTGCCATATCAACACATTCATCGAGATTCCTGTTTGCGAGCCCTGGCGTCCATTTTTTGCTCATAA
 ACGTGGAAGTCAAGATGATCAGAGCATCCTGAGTATCTGTTCTGTGATCAATCCATCCACTGCCAGTGCC
 CCAGTTTCTTCTGAGGGAAGTAA

Restriction Sites:	SgfI-MluI
ACCN:	NM_001277397
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).


[View online »](#)

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001277397.1, NP_001264326.1</u>
RefSeq Size:	864 bp
RefSeq ORF:	864 bp
Locus ID:	729974
UniProt ID:	<u>F8VTS6</u>
Cytogenetics:	19q13.42