

## Product datasheet for **SC335111**

### LETM2 (NM\_001286787) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	LETM2 (NM_001286787) Human Untagged Clone
Tag:	Tag Free
Symbol:	LETM2
Synonyms:	SLC55A2
Vector:	<u>pCMV6 series</u>
Fully Sequenced ORF:	>NCBI ORF sequence for NM_001286787, the custom clone sequence may differ by one or more nucleotides

```
ATGTTGCCATCAACTTTTGAAAGTGAATCCAAAAAGGAAGAAAAACAGAAAAAGAAAATGGCTGTAAAGT
TGGAACTAGCAAATTTCTTCAAGAAACCATGACAGAAATGGCAAGGAGGAACAGAGCCAAGATGGGCGA
TGCCCTACACAGCTCTCATCTACGTGAAGCAGGTCCAGACAGGCCACAAGCCCAGCACAAAGGAGATA
GTTTCGCTTCTCCAACTATTTGAGGACCAGCTGGCCCTGGAACACTTAGATCGCCCTCAGCTGGTTGCC
TTTGCAAAGTGTGGAATTGCAGACATTTGGAACCAACAACCTGCTCCGCTTTCAGCTCCTGATGAACT
GAAGTCTATAAAAGCAGATGATGAAATAATTGCCAAGGAAGGGGTGACAGCATTGAGTGTATCAGAACTA
CAGGCTGCCTGTAGGGCCCGAGGGATGAGATCACTGGGTCTCACGGAGGAACAACCTGCGACAACAGCTCA
CGGAGTGGCAGGACCTCCACCTGAAGGAGAACGTCCCTCCTTCCCTTTTGCTCCTGTCCCGCACCTTCTA
CCTGATAGATGTGAAGCCCAAGCCGATTGAGATACCACTCAGTGGGGAGGCTCCAAGACTGATATTCTT
GTGGAATTACCTACTTTCACTGAATCTAAAGAGAACATGGTGGATCTTGCACCTCAACTGAAGGGAAC
AGGATGAAGACTTTATACAGCCGCCACCAGTTACATCATCACCATAACACCATCAACACCTATTTTCATT
ACCTAAAGGACCCATCACTTCTTCTGAAGAACCCTACACTCCAGGCCAAATCACAATGACGGCCCGAAG
AGCAAGGCTAGTTCAAAAGGAGCATAA
```

Restriction Sites:	Sgfl-Mlul
ACCN:	NM_001286787
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).



[View online »](#)

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_001286787.1](#), [NP\\_001273716.1](#)

**RefSeq Size:** 2817 bp

**RefSeq ORF:** 867 bp

**Locus ID:** 137994

**Cytogenetics:** 8p11.23