

Product datasheet for SC335197

Mu Opioid Receptor (OPRM1) (NM_001285526) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Mu Opioid Receptor (OPRM1) (NM_001285526) Human Untagged Clone
Tag:	Tag Free
Symbol:	OPRM1
Synonyms:	LMOR; M-OR-1; MOP; MOR; MOR1; OPRM
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC335197 representing NM_001285526. Blue=Insert sequence Red=Cloning site Green=Tag(s)


```
GCTCGTTT TAGTGAACCGTCAGAATTTTGT AATACGACTCACTATAGGGCGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCC GCGATCGCC
ATGAAGACTGCCACCAACATCTACATTTTCAACCTTGCTCTGGCAGATGCCTTAGCCACCAGTACCC
CCCTTCCAGAGTGTGAATTACCTAATGGGAACATGGCCATTTGGAACCATCCTTTGCAAGATAGTGATC
TCCATAGATTACTATAACATGTTCAACAGCATATTCACCCTCTGCACCATGAGTGTGATCGATACATT
GCAGTCTGCCACCCTGTCAAGGCCTTAGATTTCCGTACTCCCGAAATGCCAAAATTATCAATGTCTGC
AACTGGATCCTCTCTTCAGCCATTGGTCTTCTGTAAATGTTTATGGCTACAACAAAATACAGGCAAGGT
TCCATAGATTGTACACTAACATTCTCTCATCCAACCTGGTACTGGGAAAACCTGCTGAAGATCTGTGTT
TTCATCTTCGCCTTATTATGCCAGTGCTCATATTACCGTGTGCTATGGACTGATGATCTTGGCCCTC
AAGAGTGTCCGCATGCTCTCTGGCTCCAAGAAAAGGACAGGAATCTTGAAGGATCACCAGGATGGTG
CTGGTGGTGGTGGCTGTGTTTCATCGTCTGCTGGACTCCCATTACATTTACGTCATCATTAAAGCCTTG
GTTACAATCCCAGAACTACGTTCCAGACTGTTTCTTGGCACTTCTGCATTGCTCTAGGTTACACAAAC
AGCTGCCTCAACCCAGTCTTTATGCATTTCTGGATGAAAACCTCAAACGATGCTTCAGAGAGTTCTGT
ATCCCAACCTCTCCAACATTGAGCAACAAAACCTCACTCGAATTCGTCAGAACACTAGAGACCACCC
TCCACGGCCAATACAGTGGATAGAACTAATCATCAGCTAGAAAATCTGGAAGCAGAACTGCTCCGTTG
CCCTAA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: Sgfl-MluI



[View online »](#)

Plasmid Map:



ACCN: NM_001285526

Insert Size: 903 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001285526.1](#)

RefSeq Size: 15825 bp

RefSeq ORF: 903 bp

Locus ID: 4988

UniProt ID: [P35372](#)

Cytogenetics: 6q25.2

Protein Families: Druggable Genome, GPCR, Transmembrane

Protein Pathways: Neuroactive ligand-receptor interaction

MW: 34.4 kDa

Gene Summary: This gene encodes one of at least three opioid receptors in humans; the mu opioid receptor (MOR). The MOR is the principal target of endogenous opioid peptides and opioid analgesic agents such as beta-endorphin and enkephalins. The MOR also has an important role in dependence to other drugs of abuse, such as nicotine, cocaine, and alcohol via its modulation of the dopamine system. The NM_001008503.2:c.118A>G allele has been associated with opioid and alcohol addiction and variations in pain sensitivity but evidence for it having a causal role is conflicting. Multiple transcript variants encoding different isoforms have been found for this gene. Though the canonical MOR belongs to the superfamily of 7-transmembrane-spanning G-protein-coupled receptors some isoforms of this gene have only 6 transmembrane domains. [provided by RefSeq, Oct 2013]

Transcript Variant: This variant (MOR-1K2) represents use of an alternate promoter and 5' UTR and uses a downstream start codon, compared to variant MOR-1i. The resulting isoform (MOR-1G1) has a shorter N-terminus, compared to isoform MOR-1i. Variants MOR-1G1, MOR-1K1, and MOR-1K2 encode the same isoform (MOR-1G1). Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.