

Product datasheet for SC335220

OR8G2P (NM_001291438) Human Untagged Clone

Product data:


Product Type:	Expression Plasmids
Product Name:	OR8G2P (NM_001291438) Human Untagged Clone
Tag:	Tag Free
Symbol:	OR8G2P
Synonyms:	HSTPCR120; OR8G2; OR8G4; ORL206; ORL486; TPCR120
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC335220 representing NM_001291438. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTT TAGTGAACCGTCAGAATTTTGT AATACGACTCACTATAGGGCGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCC CGCATCGCC
ATGGTTTTCTTCTCCGTAGAAACTGACCAAAGGAAAATGTCAGCAGGAAACCATTCTCAGTGACT
GAGTTCATTCTGGCTGGGCTCTCAGAACAGCCAGAGCTCCAGCTGCGCCTCTTCTCCTGTTCTTAGGA
ATCTATGTGGTCACAGTGGTGGCAACTGAGCATGATCACACTGATTGGGCTCAGTTCTCACCTGCAT
ACCCCATGTACTATTTCTCAGTGGTCTGTCTTCATTGATCTCTGCCATTCCACTATCATTACCCCC
AAAATGCTGGTGAACCTTGTGACAGAGAAGAACATCATCTCTACCTGAATGCATGACTCAGCTTAC
TTCTTCCTCATTTTTGCTATTGCAGAGTGCACATGTTGGCTGTAAACGGCATATGACCGCTATGTTGCC
ATCTGCAGCCCCTTGCTGTACAATGTCATCATGTCCTATCACCCTGCTTCTGGCTCACAGTGGGAGTT
TACGTTTTAGGCATCCTTGGATCTACAATTCACACCGGCTTTATGTTGAGACTCTTTTTGTGCAAGACT
AATGTGATTAACCATTATTTTTGTGATCTCTTCCCTCTCTTGGGGCTCTCCTGCTCCAGCACCTACATC
AATGAATTACTGGTCTGGTCTTGGAGTGCATTTAACATCTGACGCTGCCTTAACCATCCTTGCTTCT
TACATCTTTATCATTGCCAGCATCTCCGATTCGCTCCACTGAGGGCAGGTCAAAGCCTTCAGCACT
TGCAGCTCCACATCTTGGCTGTTGCTGTTTTCTTGGGCTGCAGCATTATGTACCTGCAGCCATCA
TCTGTAGCTCCATGGACCAGAGGAAAGTCTCTGTGTTTTATACTACTATTGTGCCCATGCTGAAC
CCCCAATCTATAGCCTAA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: SgfI-MluI



[View online »](#)

Plasmid Map:


ACCN: NM_001291438

Insert Size: 915 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001291438.1](#)

RefSeq Size: 1025 bp

RefSeq ORF: 915 bp

Locus ID: 26492

UniProt ID: [Q61F36](#)

Cytogenetics: 11q24.2

Protein Families:	Druggable Genome
Protein Pathways:	Olfactory transduction
MW:	34 kDa
Gene Summary:	<p>Olfactory receptors interact with odorant molecules in the nose, to initiate a neuronal response that triggers the perception of a smell. The olfactory receptor proteins are members of a large family of G-protein-coupled receptors (GPCR) arising from single coding-exon genes. Olfactory receptors share a 7-transmembrane domain structure with many neurotransmitter and hormone receptors and are responsible for the recognition and G protein-mediated transduction of odorant signals. The olfactory receptor gene family is the largest in the genome. The nomenclature assigned to the olfactory receptor genes and proteins for this organism is independent of other organisms. [provided by RefSeq, Jul 2008]</p>