

Product datasheet for **SC335287**

B3GAT3 (NM_001288723) Human Untagged Clone

Product data:


Product Type:	Expression Plasmids
Product Name:	B3GAT3 (NM_001288723) Human Untagged Clone
Tag:	Tag Free
Symbol:	B3GAT3
Synonyms:	GLCATI; glcUAT-I; JDSCD
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC335287 representing NM_001288723. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTGTAGTGAACCGTCAGAATTTGTAAATACGACTACTATAGGGCGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGAAGCTGAAGCTGAAGAACGTGTTTCTCGCCTACTTCTGGTGTGATCGCCGGCCTCTCTACGGC
CTGGTACAGCTCGCCAGCCATGTGACTGCCTTCTCCCTGCGGGCAGCAGCCGAGCAGCTACGGCAG
AAGGATCTGAGGATTTCCAGCTGCAAGCGAACTCCGACGGCCACCCCTGCCCTGCCAGCCCCCT
GAACCCGAGGCCCTGCCTACTATCTATGTTGTTACCCCCACCTATGCCAGGCTGGTACAGAAGGCAGAG
CTGGTACGACTGTCCAGACACTGAGCCTGGTGGCCCGGCTGCATTGGCTGCTGGTGGAGGATGCTGAG
GGTCCCACCCCGTGGTCTCAGGGCTGTGGCTGCCTCTGGCCTCTTACACACCTGGTGGTCTC
ACGCCAAAGCCAGCGGCTTCGGGAGGGCAGCCTGGCTGGGTTTCATCCCGTGGTGTGAGCAGCGG
AAACAAGGCCCTGGACTGGCTCCGGGGCAGAGGGGGTGTGTGGTGGGGAGAAGGACCCACCACCA
GGGACCAAGGAGTCGTCTACTTTGCTGACGATGACAACCTACAGCCGGGAGCTGTTTGAGGAGATG
CGCTGGACCCGTGGTGTCTCAGTGTGGCCTGTGGGGTGGTGGGGCCCTGCGATTTCGAGGGCCCTCAG
GTACAGGACGGCCGGGTAGTGGGCTTCCACACAGCATGGGAGCCAGCAGGCCCTTCCCTGTGGATAG
GCTGGATTTGCCGTGGCCCTGCCCTTGTGTTAGATAAGCCCAATGCCCAATTTGATTCCACCGCTCCC
CGGGGCCACCTGGAGAGCAGTCTTCTGAGCCACCTTGTGGATCCCAAGGACCTGGAGCCACGGGCTGCC
AACTGCACCTCGGACAGAGTCTCGTGTGTACCCAGGCTGGAGTGCAGTAG
ACGCGTACGCGGCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: SgfI-MluI



[View online »](#)

Plasmid Map:


ACCN: NM_001288723

Insert Size: 948 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001288723.1](#)

RefSeq Size: 2098 bp

RefSeq ORF: 948 bp

Locus ID: 26229

UniProt ID: [O94766](#)

Cytogenetics: 11q12.3

Protein Families:	Transmembrane
Protein Pathways:	Chondroitin sulfate biosynthesis, Heparan sulfate biosynthesis, Metabolic pathways
MW:	34.6 kDa
Gene Summary:	<p>The protein encoded by this gene is a member of the glucuronyltransferase gene family, enzymes that exhibit strict acceptor specificity, recognizing nonreducing terminal sugars and their anomeric linkages. This gene product catalyzes the formation of the glycosaminoglycan-protein linkage by way of a glucuronyl transfer reaction in the final step of the biosynthesis of the linkage region of proteoglycans. A pseudogene of this gene has been identified on chromosome 3. [provided by RefSeq, Dec 2013]</p> <p>Transcript Variant: This variant (4) uses an alternate splice site in its 3' coding region, compared to variant 1. The encoded isoform (4) has a shorter and distinct C-terminus, compared to isoform 1.</p>