

Product datasheet for **SC335317**

GMPR2 (NM_001283023) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	GMPR2 (NM_001283023) Human Untagged Clone
Tag:	Tag Free
Symbol:	GMPR2
Synonyms:	GMPR 2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC335317 representing NM_001283023. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTT TAGTGAACCGTCAGAATTTTGT AATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCC GCGATCGCC
ATGCCTCATATTGACAACGATGTGAACTGGACTTCAAGGATGTCCTTTTGAGGCCAAACGCAGTACC
CTTAAGTCTCGAAGTGAGGTGGATCTCACAAGATCCTTTTCATTCGGAACCTAAAGCAGACATACTCT
GGGGTCCCATCTTCTGCTGCCAATATGGATACTGTGGCACCTTTGAGATGGCCAAGTTCTCTGTAAG
CATCTGGCTGCCAGCTCAGGCACAGGCTCTTCTGACTTTGAGCAGCTGGAACAGATCCTGGAAGCTATT
CCCCAGGTGAAGTATATATGCCTGGATGTGGCAAATGGCTACTCTGAACACTTTGTTGAATTTGAAAA
GATGTACGGAAGCGCTTCCCCAGCACACCATCATGGCAGGGAATGGTAAACAGGAGAGATGGTAGAA
GAGCTCATCCTTTCTGGGCTGACATCATCAAAGTGGGAATTGGCCAGGCTCTGTGTACTACTCGG
AAGAAAACCTGGAGTGGGGTATCCACAGCTCAGCGCAGTGATGGAGTGTGCAGATGCTGCTCATGGCTC
AAAGGCCACATCATTTAGATGGAGTTGCAGCTGTCTGGGGATGTGGCCAAGGCTTTTGGGGCAGGA
GCTGACTTCGTGATGCTGGGTGGCATGCTGGTGGGCACAGTGAGTCAGGTGGTGAAGCTCATCGAGAGG
GATGGCAAGAAGTACAAGCTCTTCTATGGAATGAGTTCTGAAATGGCCATGAAGAAGTATGCTGGGGC
GTGGCTGAGTACAGAGCCTCAGAGGAAAGACAGTGAAGTTCCTTTAAAGGAGATGTGGAACATACC
ATCCGAGACATCCTAGGAGGGATCCGCTCTACGTGTACCTATGTGGGAGCAGCTAAGCTCAAAGAGTTG
AGCAGGAGAACTACCTTCATCCGAGTACCCAGCAGTGAATCCAATCTTCAGTGAGCGCTGC TAG
ACGCGTACGCGGCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites:	Sgfl-Mlul
ACCN:	NM_001283023
Insert Size:	963 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001283023.1
RefSeq Size:	1585 bp
RefSeq ORF:	963 bp
Locus ID:	51292
UniProt ID:	Q9P2T1
Cytogenetics:	14q12
Protein Families:	Druggable Genome
Protein Pathways:	Purine metabolism
MW:	34.6 kDa
Gene Summary:	<p>This gene encodes an enzyme that catalyzes the irreversible and NADPH-dependent reductive deamination of guanosine monophosphate (GMP) to inosine monophosphate (IMP). The protein also functions in the re-utilization of free intracellular bases and purine nucleosides. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Apr 2017]</p> <p>Transcript Variant: This variant (7) differs in the 5' UTR and contains multiple differences in the coding region, compared to variant 1, including initiation of translation at a downstream in-frame start codon. The encoded isoform (5) is shorter than isoform 1.</p>