

## Product datasheet for **SC335330**

### MYADM (NM\_001290191) Human Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** MYADM (NM\_001290191) Human Untagged Clone  
**Tag:** Tag Free  
**Symbol:** MYADM  
**Synonyms:** SB135  
**Vector:** pCMV6-Entry (PS100001)  
**Fully Sequenced ORF:** >SC335330 representing NM\_001290191.  
Blue=Insert sequence Red=Cloning site Green=Tag(s)


```
ATGCCAGTGACGGTAACCCGCACCACCATCACAACCACCACGACGTCATCTTCGGGCCTGGGGTCCCC
ATGATCGTGGGGTCCCCTCGGGCCCTGACACAGCCCCTGGGTCTCCTTCGCCTGCTGCAGCTGGTGTCT
ACCTGCGTGGCCTTCTCGCTGGTGGCTAGCGTGGGCGCCTGGACGGGTCCATGGGCAACTGGTCCATG
TTCACCTGGTGTCTTCTGCTTCTCCGTGACCCTGATCATCCTCATCGTGGAGCTGTGCGGGCTCCAGGCC
CGCTTCCCCTGTCTTGGCGCAACTTCCCATCACCTTCGCCTGCTATGCGGCCCTTCTGCTCTCG
GCCTCCATCATCTACCCACCACCTATGTCCAGTTCCTGTCCCACGGCCGTTTCGCGGGACCACGCCATC
GCCGCCACCTTCTTCTCCTGCATCGCGTGTGTGGCTTACGCCACCGAAGTGGCCTGGACCCGGGCCCGG
CCCGGCGAGATCACTGGCTATATGGCCACCGTACCCGGGCTGCTGAAGGTGCTGGAGACCTTCGTTGCC
TGCATCATCTTCGCGTTCATCAGCGACCCCAACCTGTACCAGCACCAGCCGGCCCTGGAGTGGTGCCTG
GCGGTGTACGCCATCTGCTTTCATCCTAGCGGCCATCGCCATCCTGCTGAACCTGGGGGAGTGCACCAAC
GTGCTACCCATCCCCTTCCCAGCTTCTGTGCGGGCTGGCCTTGTCTGTCTCCTCTATGCCACC
GCCCTTGTCTCTGGCCCTCTACCAGTTCGATGAGAAGTATGGCGGCCAGCCTCGGCGCTCGAGAGAT
GTAAGCTGCAGCCGACCCATGCCTACTACGTGTGTGCCTGGGACCGCCACTGGCTGTGGCCATCCTG
ACGGCCATCAACCTACTGGCGTATGTGGCTGACCTGGTGCCTCTGCCACCTGGTTTTTGTCAAGGTC
TAA
```

**Restriction Sites:** Sgfl-Mlul



[View online »](#)

## Plasmid Map:



**ACCN:** NM\_001290191

**Insert Size:** 969 bp

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_001290191.1](#)

RefSeq Size: 3083 bp  
RefSeq ORF: 969 bp  
Locus ID: 91663  
UniProt ID: [Q96S97](#)  
Cytogenetics: 19q13.42  
Protein Families: Transmembrane  
MW: 35.3 kDa