

Product datasheet for **SC335398**

SLC35A2 (NM_001282649) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	SLC35A2 (NM_001282649) Human Untagged Clone
Tag:	Tag Free
Symbol:	SLC35A2
Synonyms:	CDG2M; CDGX; UDP-Gal-Tr; UGALT; UGAT; UGT; UGT1; UGT2; UGTL
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC335398 representing NM_001282649. Blue=Insert sequence Red=Cloning site Green=Tag(s)


```
GCTCGTTT TAGTGAACCGTCAGAATTTTGT AATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCC GCGATCGCC
ATGGCAGCGGTTGGGGCTGGTGGTTCCACCGCGCGCCCGGGCCAGGGGCGGTTCCGCGGGTGCATTG
GAGCCGGGACCGCCAGTGC GGTAACGTGAAGCACCTGGTTCTCTTCTCCATGAGGCTGTCTGGTG
CAGTATGTGGACAGCTCAAGCTCGCAGTGCCCTCTCATCTACACCTGCAGAATAACCTCCAGTAT
GTTGCCATCTTAACCTACCAGTGC CACTTCCAGGTGACATACCAGCTGAAGATCCTGACCACAGCG
CTGTTCCGTGCTCATGCTGAATCGCAGCCTTCCCGGCTGCAGTGGGCTCCCTGCTGCTCCTCTT
ACTGGCGTCGCCATTGTCCAGGCACAGCAAGCCGGTGGGGAGGCCACGGCCACTGGATCAGAACCCT
GGGCAGGCCTGGCAGCGCTGTGGCCTCCTGTCTCCTCCGGCTTCGCAGGTGTCTACTTTGAGAAG
ATCCTCAAAGGCAGCTCAGGCTCCGTGTGGCTGCGCAACCTGCAACTGGGCTCTTCGGCACAGCACTG
GGCCTGGTGGGGCTCTGGTGGGCTGAGGGTACCGCCGTGGCCACCCGTGGTTTCTTTTTGGGTACACA
CCTGCTGTCTGGGGCTGGTGTCAACCAGGCCTTCGGCGGGCTACTGGTGGCTGTGGTTGTCAAGTAC
GCTGACAATATCCTCAAGGGCTTTGCCACCTCCCTGTCCATTGTGCTGTCCACTGTTGCCTCCATTCCG
CTCTTTGGCTTCCAGTGGACCCATTATTTGCCCTTGGCGCTGGACTCGTCATTGGTGTGTCTACCTC
TACAGCCTTCCCGAGGTGCAGCCAAAGCCATAGCCTCTGCCTCTGCCTCCGCCCTCCGGGCCCTGCGTT
CACCAGCAGCTCCCGGCAGCCACCACCAGCAGTGTCTCCACCCTGGAGACCTCATCACGGAG
CCTTTTCTGCCAAAGTCAGTGTGGTGAAGTGA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: SgfI-MluI



[View online »](#)

Plasmid Map:



ACCN: NM_001282649

Insert Size: 999 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001282649.1](#)

RefSeq Size: 2436 bp

RefSeq ORF: 999 bp

Locus ID: 7355

UniProt ID: [P78381](#)

Cytogenetics: Xp11.23

Protein Families: Transmembrane

MW: 34.2 kDa

Gene Summary: This gene encodes a member of the nucleotide-sugar transporter family. The encoded protein is a multi-pass membrane protein. It transports UDP-galactose from the cytosol into Golgi vesicles, where it serves as a glycosyl donor for the generation of glycans. Mutations in this gene cause congenital disorder of glycosylation type II_m (CDG2M). Multiple alternatively spliced transcript variants encoding distinct isoforms have been found for this gene. [provided by RefSeq, Oct 2014]
Transcript Variant: This variant (6) lacks an alternate in-frame exon in the 5' coding region and differs in the 3' exon structure, compared to variant 1. The encoded isoform (f) is shorter and has a distinct C-terminus, compared to isoform a.