

Product datasheet for **SC335434**

CCDC198 (NM_001283056) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	CCDC198 (NM_001283056) Human Untagged Clone
Tag:	Tag Free
Symbol:	CCDC198
Synonyms:	C14orf105
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC335434 representing NM_001283056. Blue=Insert sequence Red=Cloning site Green=Tag(s)


```
GCTCGTTTGTAGTGAACCGTCAGAATTTTGTAAACGACTACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGGCTGAGTCACTCTAAGACTCACCTTAGGGTGATCAAAGTAGCACCTTTGCAAAACAAAGAGGTA
GAGACTCCCTCGGCTGGCCGTGTGGACTTTGCATTCAATCAGAATTTGGAAGAAAAGACTTCATATTCA
CTGGCAAGACTGCAGGACCAGAATAAAGCCTTGAAGGGCAGCTGCCACCTTTACAAGAAAAGTGGTAT
GGAAGATATTCTACAGCATCCAGAGACATGTATTTTGACATCCCCTGGAACACAGAGAAACAAGTATT
ATAAAAGGCATCCACCCAAAGACTTCAAAGCTTGAACCCATTGACTTGCCACGAGTAATTACTTCA
GGAAGACTCCTGAGCCAGCGAGAAGCCAGGACAATGCACAAAGCAAAGGTAAGAAAAGAAAATGCAA
ACTCCAATGTATACTTCTGAGAACAGACAATATTTGCATAAGATGCAAGTCTGGAAATGATCCGTA
AGACAAGAGGCCCAATGGAGTTAAAGAAAAGTCTTCATGGAGAGGCAAGAATTAATAAGCAAAGTCCA
AGGGACCAATAAGCCAAGAAAACCTTCAAAGCACCCCAAGGAATGATGACCATGACCTTCTAACCATG
TTGCTGATGAAATCTTGAACAGAGGTCACGACAGAGTTTCGCTCCTGTTACCCAGGCTGGAGTGAAT
GGCATGATCTCGGCTCACCGAACTTCGCCTCCCAAGTTCAAATGATTCTCCTGCCTCAGCCTCCCGA
GTAGCTGGGATTACAGGAAATCAAAGAATACAGAATTTTGAACATCAAGCAGTGAATAACTACTGT
CCCTGGAAAATGGCAAAATGAAACATGGCTTCATGAACAAGAGGCCAGGGACAGCTTCTCTGGGAC
AGTTCAGCTCTGACTCAGATGAGCAGGGGAAAGATGAGAAGAACCCACGAGCACTGGTGAGGACCAGG
ACAGAGAGAATCCCACTTTTCGATGAGTTTTTGTCAAGAAATAA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: SgfI-MluI



[View online »](#)

Plasmid Map:



ACCN: NM_001283056

Insert Size: 1011 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001283056.1](#)

RefSeq Size: 3214 bp

RefSeq ORF: 1011 bp

Locus ID: 55195

UniProt ID: [Q9NVL8](#)

Cytogenetics: 14q22.3

MW: 38.9 kDa