

Product datasheet for **SC335491**

PPP1R7 (NM_001282412) Human Untagged Clone

Product data:


Product Type:	Expression Plasmids
Product Name:	PPP1R7 (NM_001282412) Human Untagged Clone
Tag:	Tag Free
Symbol:	PPP1R7
Synonyms:	SDS22
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC335491 representing NM_001282412. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTGTAGTAACCGTCAGAATTTTGTAAACGACTACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGTTGACAGGCGGGTCGAGTCTGAAGAATCCGGCGATGAAGAAGGGAAGAAACACAGCAGTGGCATC
GTGGCCGACCTCAGTGAACAGAGCCTGAAGGATGGGGAGGAGCGGGGGAGGAGGACCCAGAAGAAGAA
CATGAGCTGCCTGTGGACATGAAACCATCAACCTGGACAGAGATGCAGAGGATGTTGATTTGAATCAC
TATCGCATAGGGAAGATTGAAGGATTTGAGGACTGAAGAAAGTGAAGACTCTCGCTCCGCCAAAAT
TTAATTAATGCATTGAGAATCTGGAGGAGCTACAGAGTCTCGAGAGCTGGATCTTACGACAACCAG
ATCAAGAAGATTGAGAATCTGGAGGCGCTAACAGAGCTGGAGATTCTAGATATTTCTTTAATCTGCTG
AGAAACATCGAAGGGGTTGACAAGTTGACACGACTGAAAAAACTCTTCTGGTCAACAATAAAATCAGT
AAAATTGAGAACTTAAGCAACTTACATCAACTACAGATGCTAGAGCTGGGATCTAACCAGCATCCGGGCA
ATCGAAAAATATCGACACCTTAACCAACCTGGAGAGTTTGTGTTTGGGGAAAAACAAAATTAATAAATT
CAGAACCTGGATGCGCTCACCAACCTGACAGTCTCAGTATGCAGAGCAACCGGCTGACCAAGATCGAG
GGTCTGCAGAACCTGGTGAACCTGCGGGAGCTGTACCTTAGCCACAATGGCATCGAGGTCATCGAGGGC
CTGGAGAACAAATAACAACTCACGATGTTGGACATTGCATCAAAATAGAATCAAAAAGATTGAAAAATC
AGCCATTAACAGAGCTGCAAGAGTTCTGGATGAACGACAATCTCCTTGAGAGCTGGAGCGACCTCGAC
GAGCTGAAGGGAGCCAGGAGCCTGGAGACAGTGTACCTGGAGCGGAACCCCTTGCAGAAGGACCCCCAG
TACCGGGCGGAAGGTCATGCTCGCCCTCCCCTCCGTGCGGCAGATCGATGCCACGTTTCGTCAGGTTCTGA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: SgfI-MluI



[View online »](#)

Plasmid Map:


ACCN: NM_001282412

Insert Size: 1035 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001282412.1](#)

RefSeq Size: 2007 bp

RefSeq ORF: 1035 bp

Locus ID: 5510

UniProt ID: [Q15435](#)

Cytogenetics: 2q37.3

Protein Families: Druggable Genome, Phosphatase

MW: 39.8 kDa

Gene Summary: This gene encodes a protein subunit that regulates the activity of the serine/threonine phosphatase, protein phosphatase-1. The encoded protein is required for completion of the mitotic cycle and for targeting protein phosphatase-1 to mitotic kinetochores. Alternate splicing results in multiple transcript variants. [provided by RefSeq, Sep 2013]
Transcript Variant: This variant (5) uses an alternate 5' structure, lacks a portion of the 5' coding region, and initiates translation at an alternate start codon compared to variant 1. The encoded protein (isoform 5) is shorter and has a distinct N-terminus compared to isoform 1.
Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.