

Product datasheet for SC335507

ETNK2 (NM_001297762) Human Untagged Clone

Product data:


Product Type:	Expression Plasmids
Product Name:	ETNK2 (NM_001297762) Human Untagged Clone
Tag:	Tag Free
Symbol:	ETNK2
Synonyms:	EKI2; HMFT1716
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC335507 representing NM_001297762. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTT TAGTGAACCGTCAGAATTTTGT AATACGACTCACTATAGGGCGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCC GCGATCGCC
ATGGCTGTGCCCCCTTCGGCCCTCAGCCGCGCGCTCCTTTACCTGAGGAGGCACACGCCCTTGCCCG
CAGTGCTCATGGGGCATGGAGGAGAAGGCGGGCCAGCGCCAGCTGCCGGGAGCCGGGGCCCCCG
AGGGCCGCGCCGTCGCGTACTTCGGCATTTCGTGGACCCGGACGACATCCTTCCCGGGGCCCTGCGC
CTCATCCAGGAGCTGCGGCCGATTGGAAACCCGAGCAAGTTCGGACCAAGCGTTCACGGATGGCATC
ACCAACAAGCTGGTGGCCTGCTATGTGGAGGAGGACATGCAGGACTGCGTGTGGTCCGGTGTATGGG
GAGCGGACGGAGCTGCTGGTGGACCGGGAGAATGAGGTGAGAACTTCCAGCTGCTGCGAGCACACAGC
TGTGCCCCAAACTCTACTGCACCTCCAGAATGGGCTGTGCTATGAGTACATGCAGGGTGTGGCCCTG
GAGCCTGAGCACATCCGTGAGCCCGGCTTTTCAGCCTTCTGCAGATGTCCCTAAGGTAGAGGTGTTG
GAACGGGAGCTGGCTGGCTGAAGGAGCATCTGTCCCAGCTGGAGTCCCCTGTGGTGTTTTGTACAAT
GACCTGCTCTGAAGAATATCATCTATGACAGCATCAAAGGTCACGTGCGGTTTATTGACTATGAATAT
GCTGGCTACAACCTACCAAGCTTTTGACATTGGCAACCATTTCAATGAGTTTGAGGCGTGAATGAGGTG
GATTACTGCCTGTACCCGGCGGGGAGACCCAGCTGCAGTGGCTGCACTACTACCTGCAGGCACAAAAG
GGGATGGCCGTGACCCCCAGGGAGGTGCAAAGGCTCTACGTGCAAGTCAACAAGTTTGCCTGGCGTCT
CACCTTCTGGGCTCTCTGGCCCTCATCCAGAACCAGTACTCCACCATCGACTTTGATTTCTCAGG
TACGCAGTGATCCGATTCAACCAGTACTCAAGGTGAAGCCTCAAGCGTCAGCCTTGAGATGCCAAAG
TGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: SgfI-MluI



[View online »](#)

Plasmid Map:


ACCN: NM_001297762

Insert Size: 1038 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001297762.1](#)

RefSeq Size: 2515 bp

RefSeq ORF: 1038 bp

Locus ID: 55224

UniProt ID: [Q9NVF9](#)

Cytogenetics: 1q32.1

Protein Families: Druggable Genome

Protein Pathways:	Glycerophospholipid metabolism, Metabolic pathways
MW:	40 kDa
Gene Summary:	<p>The protein encoded by this gene is a member of choline/ethanolamine kinase family which catalyzes the first step of phosphatidylethanolamine (PtdEtn) biosynthesis via the cytidine diphosphate (CDP) ethanolamine pathway. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jul 2014]</p> <p>Transcript Variant: This variant (4) lacks an alternate in-frame exon, compared to variant 1. The encoded isoform (4) is shorter than isoform 1.</p>