

Product datasheet for SC335552

STK25 (NM_001282306) Human Untagged Clone


Product data:

Product Type:	Expression Plasmids
Product Name:	STK25 (NM_001282306) Human Untagged Clone
Tag:	Tag Free
Symbol:	STK25
Synonyms:	SOK1; YSK1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC335552 representing NM_001282306. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTGTAGTGAACCGTCAGAATTTTGTAAACGACTACTATAGGGCGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGAGATCGAGGACATCCAGCAGGAGATCACTGTCCCTCAAGCACCAAGCTATGGATCATCATGGAGTAC
CTGGGCGCGGCTCAGCACTGGACTTGCTAAACCAGGTCCCCTGGAGGAGACATACATTGCCACGATC
CTGCGGGAGATTCTGAAGGGCTGGATTATCTGCACTCCGAACGCAAGATCCACCGAGACATCAAAGCT
GCCAACGTGCTACTCTCGGAGCAGGGTACGCTGAAGCTGGCGGACTTTGGGGTAGCAGGGCAGCTCACA
GACACGCAGATTAAGAGGAACACATTCTGGGGCACCCCTTCTGGATGGCACCTGAGGTCATCAAGCAG
TCGGCCTACGACTTCAAGGCTGACATCTGGTCCCTGGGGATCACAGCCATCGAGCTGGCCAAGGGGGAG
CCTCCAAACTCTGACCTCCACCCATGCGCGTCTGTTCCTGATTCCCAAGAAGACGCCACCCACACTG
GAGGGCCAGCACAGCAAGCCCTTCAAGGAGTTCGTGGAGGCTGCCTCAACAAAGACCCCGATTCCGG
CCCACGGCCAAGGAGCTCCTGAAGCACAAGTTCATCACACGCTACACCAAGAAGACCTCCTTCTCACG
GAGCTCATCGACCGCTATAAGCGCTGGAAGTCAGAGGGGCATGGCGAGGAGTCCAGCTCTGAGGACTCT
GACATTGATGGCGAGGCGGAGGACGGGGAGCAGGGCCCATCTGGACGTTCCCCCTACCATCCGGCCG
AGTCCACACAGCAAGCTTCAAGGGGACGGCCCTGCACAGTTTACAGAAGCCTGCGGAGCCCGTCAAG
AGGCAGCCGAGGTCCCAGTGCCTGTCCACGCTGGTCCGGCCCGTCTTCGGAGAGCTCAAAGAGAAGCAC
AAGCAGAGCGGGGAGCGTGGGTGCGCTGGAGGAGCTGGAGAAGCCTTACAGCTGGCCGAGGAGTCC
TGCCCGGCATCTCAGACAAGCTGATGGTGCACCTGGTGGAGCGAGTGCAGAGGTTTTACACAACAGA
AACCACCTGACATCCACCCGTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: SgfI-MluI



Plasmid Map:


ACCN: NM_001282306

Insert Size: 1059 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001282306.1](#)

RefSeq Size: 2539 bp

RefSeq ORF: 1059 bp

Locus ID: 10494

Cytogenetics: 2q37.3

Protein Families: Druggable Genome, Protein Kinase

MW: 39.6 kDa

Gene Summary:

This gene encodes a member of the germinal centre kinase III (GCK III) subfamily of the sterile 20 superfamily of kinases. The encoded enzyme plays a role in serine-threonine liver kinase B1 (LKB1) signaling pathway to regulate neuronal polarization and morphology of the Golgi apparatus. The protein is translocated from the Golgi apparatus to the nucleus in response to chemical anoxia and plays a role in regulation of cell death. A pseudogene associated with this gene is located on chromosome 18. Multiple alternatively spliced transcript variants have been observed for this gene. [provided by RefSeq, Dec 2012]

Transcript Variant: This variant (3) differs in the 5' UTR and uses an alternate splice site in the 5' coding region, compared to variant 1. These differences cause translation initiation at an alternate AUG and result in an isoform (4) with a shorter and distinct N-terminus, compared to isoform 1.