

Product datasheet for SC335569

Thrombopoietin (THPO) (NM_001290028) Human Untagged Clone

Product data:


Product Type:	Expression Plasmids
Product Name:	Thrombopoietin (THPO) (NM_001290028) Human Untagged Clone
Tag:	Tag Free
Symbol:	THPO
Synonyms:	MGDF; MKCSF; ML; MPLLG; THCYT1; TPO
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC335569 representing NM_001290028. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGAGCTGACTGAATTGCTCCTCGTGGTCATGCTTCTCCTAACTGCAAGGCTAACGCTGTCCAGCCCG
GCTCCTCTGCTTGTGACCTCCGAGTCTCAGTAACTGCTTCGTGACTCCCATGTCCTTACAGCAGA
CTGAGCCAGTGCCAGAGGTTACCCCTTGGCTACACCTGCTCCTGCTGCCTGCTGTGGACTTTAGCTTG
GGAGAATGGAAAACCCAGATGGAGGAGACCAAGGCACAGGACATTCTGGGAGCAGTGACCCTTCTGCTG
GAGGGAGTGATGGCAGCACGGGGACAACCTGGGACCCACTTGCCTCTCATCCCTCCTGGGGCAGCTTCT
GGACAGGTCGGTCTCCTCCTTGGGGCCCTGCAGAGCCTCCTTGAACCCAGCTTCTCCACAGGGCAGG
ACCACAGCTACAAGGATCCCAATGCCATCTTCTGAGCTTCAAACACCTGCTCCGAGGAAAGGTGCGT
TTCTGATGCTTAGGAGGGTCCACCCTCTGCGTCAGGCGGGCCCAACCCACACAGCTGTCCCCAGC
AGAACCTCTAGTCTCACAACGAGCTCCCAAACAGGACTTCTGGATTGTTGGAGACAACTTC
ACTGCCTCAGCCAGAACTACTGGCTCTGGGCTTCTGAAGTGGCAGCAGGGATTGAGCCAAAGTTCCT
GGTCTGCTGAACCAACCTCCAGGTCCTGGACCAATCCCCGGATACCTGAACAGGATACACGAATC
TTGAATGGAACCTCGTGGACTCTTCTGGACCCCTCACGCAGGACCCCTAGGAGCCCGGACATTTCTCA
GGAACATCAGACACAGGCTCCCTGCCACCAACCTCCAGCCTGGATATTCTCTTCCCAACCCATCCT
CCTACTGGACAGTATACGCTCTTCCCTCTTCCACCCACCTTGCCACCCCTGTGGTCCAGCTCCACCCC
CTGCTTCTGACCCTTCTGCTCCAACGCCACCCCTACCAGCCCTTCTTAAACATCCTACACCCAC
TCCAGAATCTGTCTCAGGAAGGGTAA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: SgfI-MluI



[View online »](#)

Plasmid Map:


ACCN: NM_001290028

Insert Size: 1062 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001290028.1](#)

RefSeq Size: 1887 bp

RefSeq ORF: 1062 bp

Locus ID: 7066

UniProt ID: [P40225](#)

Cytogenetics: 3q27.1

Protein Families: Druggable Genome, Secreted Protein

Protein Pathways: Hematopoietic cell lineage

MW: 37.8 kDa

Gene Summary: Megakaryocytopoiesis is the cellular development process that leads to platelet production. The main functional protein encoded by this gene is a humoral growth factor that is necessary for megakaryocyte proliferation and maturation, as well as for thrombopoiesis. This protein is the ligand for MLP/C_MPL, the product of myeloproliferative leukemia virus oncogene. Mutations in this gene are the cause of thrombocythemia 1. Alternative promoter usage and differential splicing result in multiple transcript variants differing in the 5' UTR and/or coding region. Multiple AUG codons upstream of the main open reading frame (ORF) have been identified, and these upstream AUGs inhibit translation of the main ORF at different extent. [provided by RefSeq, Feb 2014]

Transcript Variant: This variant (9) represents use of the upstream promoter, but lacks an internal exon, compared to other transcripts with the same upstream promoter usage. This variant has an alternate 5' UTR exon and encodes the same isoform (1), compared to variant 1. This variant was reported in PMID: 9834204. Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.