

Product datasheet for **SC335584**

NDRG2 (NM_001282214) Human Untagged Clone

Product data:


Product Type: Expression Plasmids
Product Name: NDRG2 (NM_001282214) Human Untagged Clone
Tag: Tag Free
Symbol: NDRG2
Synonyms: SYLD
Vector: pCMV6-Entry (PS100001)
Fully Sequenced ORF: >SC335584 representing NM_001282214.
Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
ATGGCGGAGCTGCAGGAGGTGCAGATCACAGAGGAGAAGCCACTGTTGCCAGGACAGACGCCTGAGGCG
GCCAAGACTCACTCTGTGGAGACACCATACGGCTCTGTCACTTTCCTGTCTATGGCACCCCAAACCC
AAACGCCAGCGATCCTTACCTACCACGATGTGGGACTCAACTATAAATCTTGCTCCAGCCACTGTTT
CAGTTCGAGGACATGCAGGAAATCATTCAAACTTTGTGCGGGTTCATGTGGATGCCCTGGAATGGAA
GAGGGAGCCCTGTGTTCCCTTTGGGATATCAGTACCCATCTCTGGACCAGCTTGCAGACATGATCCCT
TGCCTCTGCAGTACCTAAATTTCTCTACAATAATTGGAGTTGGTGTGGAGCTGGAGCCTACATCCTG
GCGAGATATGCTCTTAACCACCCGGACACTGTTGAAGGTCTTGCTCATCAACATTGATCCCAATGCC
AAGGGTTGGATGGATTGGGCAGCCACAAGCTAACAGGCCTCACCTCTCCATTCCGGAGATGATCCTT
GGACATCTTTTCAGCCAGGAAGAGCTCTCTGAAATTCTGAGTTGATACAAAAGTACAGAAATATCATT
ACACATGCACCCAACCTGGATAACATTGAATTGTACTGGAACAGCTACAACAACCCGCCGAGACCTGAAC
TTTGAGCGTGGAGGTGATATCACCTCAGGTGCCTGTGATGCTGGTGGTAGGAGACCAAGCACCTCAT
GAAGATGCAGTGGTGAATGTAACCTAAACTGGACCCCAACCCAGACCTCGTTCCTCAAGATGGCTGAC
TCCGGAGGTGAGCCAGCTGACTCAGCCAGGCAAGCTGACCGAGGCCTCAAGTACTTCTGCAAGGC
ATGGGCTACATGGCCTCATCCTGCATGACTCGCCTGTCCCGTCTCGTACAGCCTCTGACCCAGTGCA
GCATCCGTTGATGGCAACCGGTCCGCTCTCGCACCTGTCCAGAGCAGCGAGTCTGAACTCTTTCT
TCGGGGCCCCCGGGGCACACCATGGAGGTCTCCTGTGA
```

Restriction Sites: SgfI-MluI



[View online »](#)

Plasmid Map:


ACCN: NM_001282214

Insert Size: 1074 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001282214.1](#)

RefSeq Size: 2245 bp

RefSeq ORF: 1074 bp

Locus ID: 57447

UniProt ID: [Q9UN36](#)

Cytogenetics: 14q11.2

MW: 39.3 kDa

Gene Summary: This gene is a member of the N-myc downregulated gene family which belongs to the alpha/beta hydrolase superfamily. The protein encoded by this gene is a cytoplasmic protein that may play a role in neurite outgrowth. This gene may be involved in glioblastoma carcinogenesis. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Aug 2017]

Transcript Variant: This variant (12) differs in the 5' UTR and lacks an alternate in-frame internal exon compared to variant 1. The resulting isoform (b) lacks an alternate internal segment compared to isoform a. Variants 2, 3, 5, 8, 11, and 12 all encode the same isoform (b).