

Product datasheet for **SC335641**

GDF 9 (GDF9) (NM_001288825) Human Untagged Clone

Product data:


Product Type:	Expression Plasmids
Product Name:	GDF 9 (GDF9) (NM_001288825) Human Untagged Clone
Tag:	Tag Free
Symbol:	GDF9
Synonyms:	POF14
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC335641 representing NM_001288825. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTGTAGTAACCGTCAGAATTTTGTAAATACGACTACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGAAGAAGCTCTATAAGACATATGCTACCAAGGAAGGGATTCTAAATCCAATAGAAGTCACCTCTAC
AACACTGTTTCGGCTCTTACCCCTGTACCCGGCACAAGCAGGCTCCTGGAGACCAGGTAACAGGAATC
CTTCCATCAGTGGAACTGCTATTTAACCTGGATCGCATTACTACCGTTGAACACTTACTCAAGTCAGTC
TTGCTGTACAATATCAACAACCTCAGTTTCTTTTCTCTGCTGTCAAATGTGTGCAATCTAATGATA
AAGGAGCCAAAGTCTTCTAGCAGGACTCTCGGCAGAGCTCCATACTCATTACCTTTAACTCACAGTTT
GAATTTGGAAAGAAACACAAATGGATTGAGATTGATGTGACCAGCCTCCTCAACCTTTAGTGGCCTCC
AACAAAGAGAAGTATTACATGTCTATAAATTTTACTTGCATGAAAGACCAGCTGGAGCATCCTTCAGCA
CAGAATGGTTTGTAAACATGACTCTGGTGTCCCCCTCACTGATCTTATATTTGAATGACACAAGTGCT
CAGGCTTATCACAGCTGGTATTCCCTTCACTATAAAAGGAGGCCTTCCAGGGTCTGACCAGGAGAGA
AGTCTGTCTGCCTATCCTGTGGGAGAAGAGGCTGCTGAGGATGGGAGATCTTCCCATCACCGTCACCGC
AGAGGTCAGGAACTGTCAGTTCTGAATTGAAGAAGCCCTTGGGCCAGCTTCTTCAATCTGAGTGAA
TACTTCAGACAATTTCTTCTTCCCAAATGAGTGTGAGCTCCATGACTTTAGACTTAGCTTTAGTCAG
CTGAAGTGGGACAACCTGGATTGTGGCTCCGCACAGGTACAACCTCGATACTGTAAGGGGACTGTCCA
AGGGCAGTTGGACATCGGTATGGCTCTCAGTTCACACCATGGTACAGAACATCATCTATGAGAAGCTG
GACTCCTCAGTGCCAAGACCGTCATGTGTACCTGCCAAATACAGCCCTTGAGTGTGTTTGACCATTGAG
CCCAGTGGCTCAATTGCCTATAAAGAGTACGAAGATATGATAGCTACAAAGTGCACCTGTCGTTAA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: SgfI-MluI



[View online »](#)

Plasmid Map:


ACCN: NM_001288825

Insert Size: 1101 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001288825.2](#)

RefSeq Size: 1948 bp

RefSeq ORF: 1101 bp

Locus ID: 2661

Cytogenetics: 5q31.1

Protein Families:	Adult stem cells, Cancer stem cells, Druggable Genome, Embryonic stem cells, ES Cell Differentiation/IPS, Secreted Protein, Stem cell - Pluripotency, Stem cell relevant signaling - TGFb/BMP signaling pathway, Transmembrane
MW:	41.8 kDa
Gene Summary:	<p>This gene encodes a secreted ligand of the TGF-beta (transforming growth factor-beta) superfamily of proteins. Ligands of this family bind various TGF-beta receptors leading to recruitment and activation of SMAD family transcription factors that regulate gene expression. The encoded preproprotein is proteolytically processed to generate each subunit of the disulfide-linked homodimer. This protein regulates ovarian function. Reduced expression of this gene may be associated with polycystic ovary syndrome and mutations in this gene may be more common in mothers of dizygotic twins. [provided by RefSeq, Jul 2016]</p> <p>Transcript Variant: This variant (3) contains two alternate exons in the 5' UTR, lacks a portion of the 5' coding region, and initiates translation at a downstream in-frame start codon, compared to variant 1. The encoded isoform (2) has a shorter N-terminus than isoform 1. Variants 2, 3, 4, 5, and 6 encode the same isoform (2).</p>