

Product datasheet for **SC335649**

TSNAXIP1 (NM_001288992) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	TSNAXIP1 (NM_001288992) Human Untagged Clone
Tag:	Tag Free
Symbol:	TSNAXIP1
Synonyms:	TXI1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Restriction Sites:	Sgfl-MluI
ACCN:	NM_001288992
Insert Size:	1101 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_001288992.1</u>


[View online »](#)

RefSeq Size: 2249 bp

RefSeq ORF: 1101 bp

Locus ID: 55815

UniProt ID: [Q2TAA8](#)

Cytogenetics: 16q22.1

MW: 42.3 kDa

Gene Summary: Possible role in spermatogenesis.[UniProtKB/Swiss-Prot Function]
Transcript Variant: This variant (4) differs in its 5' UTR and uses a downstream start codon, compared to variant 1. The encoded isoform (d) has a shorter N-terminus, compared to isoform a. Variants 4, 5, and 6 encode the same isoform.