

Product datasheet for **SC335663**

SEC13L1 (SEC13) (NM_001136026) Human Untagged Clone

Product data:


Product Type:	Expression Plasmids
Product Name:	SEC13L1 (SEC13) (NM_001136026) Human Untagged Clone
Tag:	Tag Free
Symbol:	SEC13
Synonyms:	D3S1231E; npp-20; SEC13L1; SEC13R
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC335663 representing NM_001136026. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGAGAGAACCTGTGCTTACTTGGTGTGTACCTCTTGAAGTCTTGCAGCCATCCTCTTCCACTGTCT
GCCTTTCTGAAATCACAAGTCAAACCTTTATACCTACAGGGCTTGTGCTGGAAAAGATGAAATGGGGAAA
ATGGTGTGAGTAATTAACACTGTGGATACCTCCCATGAGGACATGATTCACGACGCCAGATGGACTAC
TATGGCACCCGCCTGGCAACCTGCTCATCAGACAGGTCGGTCAAATCTTTGATGTGCGCAATGGAGGG
CAGATCCTTATCGCCGACCTCAGGGTCATGAGGGTCTGTGTGGCAAGTGGCCTGGGCTACCCCATG
TACGGCAACATCCTGGCATCGTGTCTCTATGACCGGAAAGTCATTATCTGGAGAGAGGAAAACGGCACC
TGGGAGAAGAGCCACGAGCATGCCGGACACGACTCCTCAGTGAATCGGTGTGCTGGGCCCCCATGAC
TACGGCCTGATCCTGGCCTGTGGGAGCTCGGATGGGGCCATCTCCCTGCTGACTTACACCGGGGAAGGC
CAATGGGAAGTAAAGAAGATCAACAACGCTCACACCATGGCTGCAATGCCGTGAGCTGGGCCCTGCT
GTTGTACTGGAAAGCCTCATAGACCACCCATCGGGGCAGAAACCAATTACATCAAGAGGTTTGCATCA
GGTGGCTGTGACAACCTCATCAAGCTGTGGAAGGAGGAGGAGGACGGCCAGTGAAGGAGGAGCAGAAG
CTAGAAGCGCACAGTACTGGGTTTCGAGATGTGGCCTGGGCCCTCCATCGGCCTGCCACCAGCACC
ATCGCCAGCTGCTCCAGGATGGTCTGTGTTTCAATTTGGACCTGTGATGATGCCTCAAGCAATACGTGG
TCCCCTAAATTGTTGCACAAGTCAACGATGTGGTGTGGCATGTGAGCTGGTCCATCACAGCCAACATC
CTGGCTGTCTCTGGTGGAGACAATAAGGTGACCCTGTGGAAGGAGTCAAGTGTGATGGGCAAGTGGGTGTGC
ATCAGTGTGATCAACAAGGGCCAGGGCTCCGTATCAGCATCAGTACAGAGGGCCAGCAGAACGAGCAG
TGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: SgfI-MluI



[View online »](#)

Plasmid Map:


ACCN: NM_001136026

Insert Size: 1107 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001136026.2](#)

RefSeq Size: 1615 bp

RefSeq ORF: 1107 bp

Locus ID: 6396

UniProt ID: [P55735](#)

Cytogenetics: 3p25.3

MW: 40.7 kDa

Gene Summary:

The protein encoded by this gene belongs to the SEC13 family of WD-repeat proteins. It is a constituent of the endoplasmic reticulum and the nuclear pore complex. It has similarity to the yeast SEC13 protein, which is required for vesicle biogenesis from endoplasmic reticulum during the transport of proteins. Multiple alternatively spliced transcript variants have been found. [provided by RefSeq, Oct 2008]

Transcript Variant: This variant (3) represents the longest transcript and encodes the longest isoform (3).