

Product datasheet for SC335675

Heme oxygenase 2 (HMOX2) (NM_001286267) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Heme oxygenase 2 (HMOX2) (NM_001286267) Human Untagged Clone
Tag:	Tag Free
Symbol:	HMOX2
Synonyms:	HO-2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC335675 representing NM_001286267. Blue=Insert sequence Red=Cloning site Green=Tag(s)


```
GCTCGTTTGTAGTAACCGTCAGAATTTTGTAAACGACTACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGAGTCTCGTCTGTTGCCAGGCTGGAGTGCAGTGGATCAGTCTAGCTCATTGCAGCCTCCACCTC
CTGGGTTCAAGAGATTCTCCTGCCTCAGCCTCCTGTGTAGCTGGGATTACAGGACCAGAGGAGCGAGAG
CAGCAAGAACCACCCAGCAGCAATGTCAGCGGAAGTGGAAACCTCAGAGGGGTAGACGAGTCAGAA
AAAAAGAACTCTGGGGCCCTAGAAAAGGAGAACCAAATGAGAATGGCTGACCTCTCGGAGCTCCTGAAG
GAAGGGACCAAGGAAGCACACGACCGGGCAGAAAACCCAGTTTGTCAAGGACTCTTGAAGGCAAC
ATTAAGAAGGAGCTGTTAAGCTGGCCACCACGGCACTTTACTTCACATACTCAGCCCTCGAGGAGGAA
ATGGAGCGCAACAAGGACCATCCAGCCTTTGCCCTTTGACTTCCCATGGAGCTGCACCGGAAGGAG
GCGCTGACCAAGGACATGGAGTATTTCTTTGGTAAAACCTGGGAGGAGCAGGTGCAGTGCCTCAAGGCT
GCCCAGAAGTACGTGGAGCGGATCCACTACATAGGGCAGAACGAGCCGGAGCTACTGGTGGCCCATGCA
TACACCCGCTACATGGGGATCTCTCGGGGGCCAGGTGCTGAAGAAGGTGGCCAGCGAGCACTGAAA
CTCCCCAGCACAGGGGAAGGGACCCAGTTCTACCTGTTTGAAGTGTGGACAATGCCCAGCAGTTCAAG
CAGCTCTACCGGGCCAGGATGAACGCCCTGGACCTGAACATGAAGACCAAGAGAGGATCGTGGAGGAG
GCCAACAAGGCTTTTGTAGTATAACATGCAGATATTCAATGAACTGGACCAGGCCGCTCCACACTGGCC
AGAGAGACCTTGGAGGATGGGTTCCCTGTACACGATGGGAAAGGAGACATGCGTAAATGCCCTTTCTAC
GCTGCTGAACAAGACAAAGGTGCCCTGGAGGGCAGCAGCTGTCCCTTCCGAACAGCTATGGCTGTGCTG
AGGAAGCCAGCCTCCAGTTCATCCTGGCCGCTGGTGTGGCCCTAGCTGCTGGACTCTTGGCCTGGTAC
TACATGTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: SgfI-MluI



[View online »](#)

Plasmid Map:



ACCN: NM_001286267

Insert Size: 1113 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001286267.1](#)

RefSeq Size: 1875 bp

RefSeq ORF: 1113 bp

Locus ID: 3163

UniProt ID: [P30519](#)

Cytogenetics: 16p13.3

Protein Families: Transmembrane

Protein Pathways: Porphyrin and chlorophyll metabolism

MW: 41.7 kDa

Gene Summary: Heme oxygenase, an essential enzyme in heme catabolism, cleaves heme to form biliverdin, which is subsequently converted to bilirubin by biliverdin reductase, and carbon monoxide, a putative neurotransmitter. Heme oxygenase activity is induced by its substrate heme and by various nonheme substances. Heme oxygenase occurs as 2 isozymes, an inducible heme oxygenase-1 and a constitutive heme oxygenase-2. HMOX1 and HMOX2 belong to the heme oxygenase family. Several alternatively spliced transcript variants encoding three different isoforms have been found for this gene. [provided by RefSeq, Oct 2013]
Transcript Variant: This variant (5) encodes the longest isoform (a).