

Product datasheet for **SC335680**

ZNF596 (NM_001287399) Human Untagged Clone

Product data:


Product Type:	Expression Plasmids
Product Name:	ZNF596 (NM_001287399) Human Untagged Clone
Tag:	Tag Free
Symbol:	ZNF596
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC335680 representing NM_001287399. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGAGAACGAAACACTTTGTAAGCAAAAAGTTGGGAAAATCTTCAGTGACTGGTTATCCTTTAATCAA
CACAAAGGAAATTCACACCAAAATGTAATTCATATGGAAGTCATCTATTTGATTATGCCTTTATCCAAAAC
TCTGCCCTTAGACCACACAGTGTGACTCACACTAGAGAGATAACATTGGAATGTCGTGTGTGTGGGAAA
ACCTTTAGCAAAAATTCTAATCTTAGGGCAGATGAGATGATTCACACTGGAGAGAAACCACCGGATGT
CATCTATGTGGGAAAGCCTTTACTCATTGCTCTGATCTTCGAAAACATGAGAGAAGTACACTGGAGAG
AAGCCATATGGATGTCATCTATGTGGGAAAGCCTTCAGTAAAAGTTCTAACCTTAGACGACATGAGATG
ATTCACACTAGAGAAAAGCACAGATATGCCATCTATGTGGGAAAGCCTTCACTCATTGCTCTGACCTT
AGAAAACATGAGAGAAGTCACTTAGGAGATAAACCATATGGATGTCTCCTATGTGGGAAAGCCTTTAGT
AAATGTTCTTACCTTAGACAACATGAAAGAAGTCACAATGGAGAGAAACCATATGAATGTCATCTATGT
GGAAAAGCCTTCTCTCATTGTTCTCACCTTAGACAACATGAGCGAAGTCACAATGGAGAGAAACCACAT
GGATGTCATCTATGTGGGAAAGCATTCACTGAATCTTCTGTGCTTAAACGACATGAGAGAATTACACT
GGAGAGAAACCATATGAGTGCCATGTATGTGGGAAAGCCTTCACTGAATCTTCTGACCTCAGACGACAT
GAGAGAAGTCACTGGAGAAAACCATATGAATGCCATCTATGCGGAAAAGCCTTCAATCACTCTTCT
GTCCTTAGACGACATGAGAGAAGTCACTGGAGAGAAACCATATGAATGCAATATATGTGGTAAAGCC
TTCAATAGAAGTTACAACCTTAGACTTCATAGAAGAGTTCACACTGGAGAGAAACCATATGTATGTCCT
CTATGTGGGAAAGCCTTTAGTAAATTTTTAACTTAGACAACATGAGAGAAGTCACTAAAAAAGCA
ATGAATATGTAA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: Sgfl-MluI



[View online »](#)

Plasmid Map:


ACCN: NM_001287399

Insert Size: 1116 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001287399.1](#)

RefSeq Size: 2779 bp

RefSeq ORF: 1116 bp

Locus ID: 169270

UniProt ID: [Q8TC21](#)

Cytogenetics: 8p23.3

Protein Families: Druggable Genome, Transcription Factors

MW: 43.3 kDa

Gene Summary: May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function]
Transcript Variant: This variant (7) differs in the 5' UTR and the 5' coding region and initiates translation at a downstream start codon, compared to variant 1. It encodes isoform 2, which is shorter at the N-terminus, compared to isoform 1. Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.