

Product datasheet for SC335706

XKR9 (NM_001287259) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	XKR9 (NM_001287259) Human Untagged Clone
Tag:	Tag Free
Symbol:	XKR9
Synonyms:	XRG9
Vector:	pCMV6-Entry (PS100001)
Fully Sequenced ORF:	>SC335706 representing NM_001287259. Blue=Insert sequence Red=Cloning site Green=Tag(s)

ATGAAATATACTAAACAGAATTTTATGATGTCAGTTCTTGGCATTATAATCTACGTAACGTGATTTAATT
 GTGGACATATGGGTATCTGTCAGATTTTCCATGAAGGACAGTATGTTTTAGTGCTTTAGCGTTAAGC
 TTTATGCTTTTGAACACTTGTGGCTCAGTGTTTAGTTATTCTTGGTTCAAGGCTGATTTAAAGAAA
 GCAGGCCAAGAAAGTCAGCATTGTTTTCTTCTACTTCATTGCTTGCAAGGAGGATTTTACAAGGTAT
 TGGTTTGCCTTAAAAAGGGTTACCATGCAGCTTTTAAATATGACAGCAATACTAGTAACCTCGTGGAA
 GAACAAATTGATCTACATAAAGAAGTTATAGATAGAGTGACTGATTTGAGCATGCTCAGACTATTTGAG
 ACCTACCTGGAAGGCTGCCCACTTATTCTTCAACTCTACATTCTTCTGGAGCATGGACAAGCGAAT
 TTCAGTCAGTATGCGGCCATCATGGTCTCTTGCTGTGCTATTTCTTGGTCAACTGTTGATTATCAAGTA
 GCTTTAAGAAAATCCTTGCTGACAAAAAGCTTCTTAATGGATTATGTCCCAAAATCACATATCTCTTT
 TACAAGTTGTTTACATTATTATCGTGGATGCTGAGTGTTGTACTTCTACTATTCTTAAATGTTAAGATT
 GCTTTATTTCTGTTGTTATTTCTTTGGTTGTTAGGTATAATATGGGCATTTAAAAACAACACCCAGTTT
 TGTACTTGATAAGTATGGAATTCCTATATAGGATTGTTGTTGGATTCACTCTTATCTTTACATTTTTT
 AATATTAAGGGACAGAATACCAAGTGCCAATGTCTGTTATTATATTGTTAGGGTACTGGGCACCTTTG
 GGGATATTGACTGTATTCTGGGTTTGGCCCTCACTATTTTAAATCCAGACTATTTTATACCTATCAGT
 ATAAGTATAGTTCTTACTCTTCTTCTTGAATTCCTTTCTTATTGTTTATTATGGGAGTTTTCACCCA
 AACAGAAGTGCAGAAACAAAATGTGATGAAATTGATGGAAAACCAGTTCTAAGAGAATGTAGAATGAGA
 TATTTCTAATGGAA**TAA**

Restriction Sites:	SgfI-MluI
ACCN:	NM_001287259
Insert Size:	1122 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).


[View online »](#)

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001287259.1</u>
RefSeq Size:	3301 bp
RefSeq ORF:	1122 bp
Locus ID:	389668
UniProt ID:	<u>Q5GH70</u>
Cytogenetics:	8q13.3
MW:	43.4 kDa