

Product datasheet for SC335731

SERPINB6 (NM_001297700) Human Untagged Clone

Product data:


Product Type:	Expression Plasmids
Product Name:	SERPINB6 (NM_001297700) Human Untagged Clone
Tag:	Tag Free
Symbol:	SERPINB6
Synonyms:	CAP; DFNB91; MSTP057; PI-6; PI6; PTI; SPI3
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC335731 representing NM_001297700. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTGTAAACGACTACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGATGTTCTCGCAGAAGCAAATGGCACCTTTGCCTTAAACCTTTGAAAACGCTGGGTAAGACAAC
TCGAAGAATGTGTTTTCTCACCATGAGCATGCTCCTGTGCCCTGGCCATGGTCTACATGGGGGCAAG
GGAAACACCGCTGCACAGATGGCCAGATACTTTCTTTCAATAAAAGTGGCGGTGGTGGAGACATCCAC
CAGGGCTTCCAGTCTTCTCACCGAAGTGAACAAGACTGGCACGCAGTACTTGCTTAGGATGGCCAAC
AGGCTCTTTGGGAAAAGTCTTGTGATTTCTCTCATCTTTAGAGATTCCTGCCAAAAATCTACCAA
GCAGAGATGGAGGAGCTTGACTTTATCAGCGCCGTAGAGAAGTCCAGAAAACACATAAACACCTGGGTA
GCTGAAAAGACAGAAGGTAAAATTGCGGAGTTGCTCTCTCCGGGCTCAGTGGATCCATTGACAAGGCTG
GTTCTGGTGAATGCTGTCTATTTTCAGAGGAAACTGGGATGAACAGTTTGACAAGGAGAACACCGAGGAG
AGACTGTTTAAAGTCAGCAAGAATGAGGAGAAACCTGTGCAAAATGATGTTTAAAGCAATCTACTTTTAAAG
AAGACCTATATAGGAGAAATATTTACCCAAATCTTGGTGCTTCCATATGTTGGCAAGGAACCTGAATATG
ATCATCATGCTTCCGGACGAGACCCTGACTTGAGAACGGTGGAGAAAGAACTCACTTACGAGAAGTTC
GTAGAAATGGACGAGGCTGGACATGATGGATGAAGAGGAGGTGGAAGTGTCCCTCCCGCGGTTTAACTA
GAGGAAAGCTACGACATGGAGAGTGTCTGCGCAACCTGGGCATGACTGATGCCTTCGAGCTGGGCAAG
GCAGACTTCTCTGGAATGTCCCAGACAGACCTGTCTGTCCAAGTCTGCACAAGTCTTTTGTGGAG
GTCAATGAGGAAGGCACGAGGCTGCAGCCGCCACAGCTGCCATCATGATGATGCGGTGTGCCAGATTC
GTCCCGGCTTCTGCGCCGACCACCCCTTCTTTTCTTCCATCCAGCACAGCAAGCAACGGGATTCTC
TTCTGCGGCCGCTTTTCTCTCCGTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: SgfI-MluI



[View online »](#)

Plasmid Map:


ACCN: NM_001297700

Insert Size: 1131 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001297700.1](#)

RefSeq Size: 1873 bp

RefSeq ORF: 1131 bp

Locus ID: 5269

UniProt ID: [P35237](#)

Cytogenetics: 6p25.2

Protein Families: Druggable Genome

MW: 42.6 kDa

Gene Summary:

The protein encoded by this gene is a member of the serpin (serine proteinase inhibitor) superfamily, and ovalbumin(ov)-serpin subfamily. It was originally discovered as a placental thrombin inhibitor. The mouse homolog was found to be expressed in the hair cells of the inner ear. Mutations in this gene are associated with nonsyndromic progressive hearing loss, suggesting that this serpin plays an important role in the inner ear in the protection against leakage of lysosomal content during stress, and that loss of this protection results in cell death and sensorineural hearing loss. Alternatively spliced transcript variants have been found for this gene. [provided by RefSeq, Sep 2010]

Transcript Variant: This variant (8) contains an alternate 5' UTR exon, compared to variant 1. Variants 1, 5, 6, 7 and 8 encode the same isoform (a).