

Product datasheet for **SC335753**

UBL7 (NM_001286740) Human Untagged Clone

Product data:


Product Type: Expression Plasmids
Product Name: UBL7 (NM_001286740) Human Untagged Clone
Tag: Tag Free
Symbol: UBL7
Synonyms: BMSC-UbP; TCBA1
Vector: pCMV6-Entry (PS100001)
Fully Sequenced ORF: >SC335753 representing NM_001286740.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

ATGTCTCTCTCAGACTGGCACCTGGCGGTGAAGCTGGCTGACCAGCCACTTACTCCAAAGTCTATTCTT
 CGGTTGCCAGAGACAGAAGTGGGAGAATACTCGCTAGGGGGCTATAGTATTTTCTGAAGCAGCTT
 ATTGCTGGCAAACCTCCAGGAGTCTGTTCCAGACCCTGAGCTGATTGATCTGATCTACTGTGGTGGGAA
 CTAAGATGACCAGACACTTGACTTCTATGGCATTCAACCTGGGTCCACTGTCCATGTTCTGCGAAAG
 TCCTGGCCTGAACCTGATCAGAAACCGGAACCTGTGGACAAAGTGGCTGCCATGAGAGAGTCCGGGTG
 TTGCACACTGCCCTGCACAGCAGCTCCTTTACAGGGAGGCGGTCTTTAAGATGCTCAGCAATAAGGAG
 TCTCTGGATCAGATCATTGTGGCCACCCAGGCCCTCAGCAGTGACCCTATTGCTCTTGGGTTCTCCAG
 GACAAGGACCTCTTCTCTGTCTTCGCTGATCCCAATATGCTTGATACGTTGGTGCCTGCTCACCCAGCC
 CTCGTCATGCCATTGTCTGGTCTGCACTCCGTAGCAGGCAGTGCCCAATGCCTGGGACTGACTCC
 TCTTCCCGGAGCATGCCCTCCAGCTCATACCGGGATATGCCAGGTGGCTTCTGTTTGAAGGGCTCTCA
 GATGATGAGGATGACTTTCACCCAAACACCAGGTCCACACCCTTAGCAGTACTCCAGCTCCCGCCA
 GCCTCCCTGGGGTACAGTGGAGCTGCTGGGCCCGGCCATCACCCAGAGTGAGCTGGCCACCGCCTTG
 GCCCTGGCCAGCACTCCGGAGAGCAGCTCTCACACCCGACTCCTGGCACCCAGGGTCAATCCTCAGGG
 ACCTACCAATGTCCTCTGGTGTCCAGTCAGGGACGCCATCACCAATGATCTCTCAGCCAAGCCCTA
 CAGCATGCCCTTCAGGCCCTCTGGGCAGCCAGCCTTCAGAGCCAGTGGCAGCCCCAGCTGCAGCAGCTA
 CGTGACATGGGCATCCAGGACGATGAGCTGAGCCTGCGGGCCCTGCAGGCCACCGGTGGGACATCCAA
 GCAGCCCTGGAGCTCATCTTTGCTGGAGGAGCCCCATGA

Restriction Sites: SgfI-MluI



[View online »](#)

Plasmid Map:


ACCN: NM_001286740

Insert Size: 1143 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001286740.1](#)

RefSeq Size: 1335 bp

RefSeq ORF: 1143 bp
Locus ID: 84993
UniProt ID: Q96S82
Cytogenetics: 15q24.1
Protein Families: Druggable Genome
MW: 40.5 kDa