

Product datasheet for SC335796

GOT2 (NM_001286220) Human Untagged Clone

Product data:


Product Type:	Expression Plasmids
Product Name:	GOT2 (NM_001286220) Human Untagged Clone
Tag:	Tag Free
Symbol:	GOT2
Synonyms:	DEE82; KAT4; KATIV; KYAT4; mitAAT
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC335796 representing NM_001286220. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTGTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGCCCTGCTGCACTCCGGCCGCGTCCTCCCGGGATCGCCGCCGCTTCCACCCGGGCCCTCGCCGCC
GCGGCCTCTGCCAGAGCCAGCTCCTGGTGGACCCATGTGGAAATGGGACCTCCAGATCCCATTCTGGGA
GTCCTGAAGCCTTAAAGAGGGACACCAATAGCAAAAAGATGAATCTGGGAGTTGGTGCCTACCGGGAT
GATAATGGAAAGCCTTACGTTCTGCCTAGCGTCCGCAAGTTTGTCACTGTGCAGACCATTCTGGAAC
GGAGCCTAAGGATCGGAGCCAGTTTTCTGCAAAGATTTTTAAGTTAGCCGAGATGTCTTTCTGCC
AAACCAACCTGGGAAACCACACCCCATCTTCAGGGATGCTGGCATGCAGCTACAAGTTATCGGTAT
TATGACCCCAAGACTTGCAGTTTGTACTTACAGGCGCTGTGGAGGATATTTCAAAAATACCAGAGCAG
AGTGTCTTCTTCTGCATGCCTGCGCCACAATCCCACGGGAGTGGACCCGCGTCCGGAACAGTGGAA
GAAATAGCAACAGTGGTGAAGAAAAGGAATCTCTTTGCGTTCTTTGACATGGCCTACCAAGGCTTTGCC
AGTGGTATGGTATAAGGATGCCTGGGCTGTGCGCCACTTATCGAACAGGGCATTATGTTTGCCTC
TGCCAATCATATGCCAAGAACATGGGCTTATATGGTGAGCGTGTAGGAGCCTTCACTATGGTCTGCAAA
GATGCGGATGAAGCCAAAAGGGTAGAGTCACAGTTGAAGATCTTGATCCGTCCTATGTTCCAAACCT
CCCCTCAATGGGGCCCGGATTGCTGCTGCCATTCTGAACACCCAGATTTGCGAAAACAATGGCTGCAA
GAAGTAAAGTCATGGCTGACCGCATCATTGGCATGCGGATCACTGGTCTCCAACCTCAAGAAGGAG
GGTTCCACCCACAATTGGCAACACATCACCGACCAAAATTGGCATGTTCTGTTTACAGGGCTAAAGCCT
GAACAGGTGGAGCGGCTGATCAAGGAGTTCTCCATCTACATGACAAAAGATGGCCGCATCTCTGTGGCA
GGGGTACCTCCAGCAACGTGGGCTACCTTGCCCATGCCATTACCAGGTACCAAGTAA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: SgfI-MluI



[View online »](#)

Plasmid Map:


ACCN: NM_001286220

Insert Size: 1164 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001286220.1](#)

RefSeq Size: 2359 bp

RefSeq ORF: 1164 bp

Locus ID: 2806

UniProt ID: [P00505](#)

Cytogenetics: 16q21

Protein Families: Stem cell - Pluripotency

Protein Pathways:	Alanine, aspartate and glutamate metabolism, Arginine and proline metabolism, Cysteine and methionine metabolism, Metabolic pathways, Phenylalanine, tyrosine and tryptophan biosynthesis, Phenylalanine metabolism, Tyrosine metabolism
MW:	43 kDa
Gene Summary:	<p>Glutamic-oxaloacetic transaminase is a pyridoxal phosphate-dependent enzyme which exists in cytoplasmic and inner-membrane mitochondrial forms, GOT1 and GOT2, respectively. GOT plays a role in amino acid metabolism and the urea and tricarboxylic acid cycles. The two enzymes are homodimeric and show close homology. Two transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Oct 2013]</p> <p>Transcript Variant: This variant (2) lacks an alternate in-frame exon compared to variant 1. The resulting isoform (2) has the same N- and C-termini but is shorter compared to isoform 1.</p>