

Product datasheet for **SC335803**

ZNF670 (NM_001204220) Human Untagged Clone

Product data:


Product Type:	Expression Plasmids
Product Name:	ZNF670 (NM_001204220) Human Untagged Clone
Tag:	Tag Free
Symbol:	ZNF670
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC335803 representing NM_001204220. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGATTCAGTGTCAATTTGAAGATGTGGCTGTGGCCTTTACTCAGGAGGAGTGGGCTTTGCTGGATCCT
TCTCAAAAGAATCTCTACAGAGATGTGATGCAAGAAATCTTCAGGAACCTGGCTTCTGTAGGAAACAAA
TCAGAAGACCAGAATATCCAAGATGACTTCAAAAATCCTGGGAGAAAATCTAAGTCATGTGGTAGAGAGA
CTGTTTGAATTAAGAAGGCAGTCAATATGGAGAAACCTTCAGCCAGGATTCAAATTTGAATCTGAAT
AAGAAAGTTTCTACTGGAGTAAAGCCATGTGAATGCAGTGTGTGGAAAAGTCTTCATATGTCATTCA
GCCCTTCATAGGCACATCCTGTCTCACATTGAAACAAACTATTTGAGTGTGAGGAATGCCAGAGAAG
TTATATCATTGCAAAACAATGTGGGAAAGCCTTTATCTCTCACCAGTGTGACAGACACATGGTAACA
CACACTAGTAATGGGCCTTATAAGGGTCCAGTGTATGAGAAGCCTTTTGATTTTCCTAGTGTATTTCAA
ATGCCTCAGAGCACTTACACTGGAGAGAAAACATATAAATGTAACATTGTGATAAAGCCTTCAATTAT
TCAAGTTATCTTCGTGAACATGAAAGAACTCATACTGGAGAGAAACCTATGCATGTAAGAAATGTGGT
AAATCATTCACTTTTCCAGTCTCTTCGCCAACATGAAAGATCTCATACTGGAGAGAAACCTATGAA
TGTAAGGAATGTGGCAAAGCCTTCAGTCGTTCCACTTACTGGGAATACATGAAAGAACGCATACTGGA
GAAAAACCCTATGAATGTATAAAATGTGGCAAAGCCTTTAGATGTTCCAGAGTCCCTCAGAGTCCATGAA
AGGACTCACAGTGGAGAAAAGCCCTATGAATGTAACAATGTGGTAAAGCCTTCAAATATTTCTAGTAAC
CTATGTGAGCATGAAAGAACTCACACTGGAGTGAACCTTATGGATGTAAGGAATGTGGTAAAGTCGTTT
ACTTCTCCAGTGCCTTCGAAGCCATGAAAGGACTCATACTGGAGAAAAACCCTATGAATGTAAGAAA
TGTGGTAAAGCCTTCAGTTGTTCCAGTTCCTTCGAAAGCATGAAAGAGCTTATATGTGGTAA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: Sgfl-MluI



[View online »](#)

Plasmid Map:


ACCN: NM_001204220

Insert Size: 1167 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001204220.1](#)

RefSeq Size: 4195 bp

RefSeq ORF: 1167 bp

Locus ID: 93474

Cytogenetics: 1q44

Protein Families: Transcription Factors

MW: 44.5 kDa

Gene Summary:

May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function]

Transcript Variant: This variant (2) uses an alternate in-frame splice site in the 3' coding region, compared to variant 1, resulting in an isoform (2) that is 1 aa shorter than isoform 1.

Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.