

Product datasheet for **SC335806**

TMPRSS5 (NM_001288752) Human Untagged Clone

Product data:


Product Type:	Expression Plasmids
Product Name:	TMPRSS5 (NM_001288752) Human Untagged Clone
Tag:	Tag Free
Symbol:	TMPRSS5
Synonyms:	SPINESIN
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC335806 representing NM_001288752. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGAGCCTGATGCTGGATGACCAACCCCTATGGAGGCCAGTATGCAGAGGAGGGCCAGGACCTGGG
ATCTTCAGAGCAGAGCCTGGAGACCAGCAGCATCCCATTCTCAGGCAGTGTGCTGGCGTTCATGCGA
CGTGGCTGTGCAGTGTGGGAGCCCTGGGGCTGCTGGCCGGTGCAGGTGTTGGCTCATGGCTCCTAGTG
CTGTATCTGTGCTGCTGCCTCTCAGCCATTCCGGGACCTTGCAGGATGAGGAGATAACTTTGAGC
TGCTCAGAGGCCAGCGCTGAGGAAGCTCTGCTCCCTGCACTTCCAAAACAGTATCTTTCAGAATAAAC
AGCGAAGACTTCTTGCTGGAAGCGCAAGTGAGGGATCAGCCACGCTGGCTCCTGGTCTGCCATGAGGGC
TGGAGCCCCGCCCTGGGGCTGCAGATCTGCTGGAGCCTTGGGCATCTCAGACTCACTCACCACAAGGGA
GTA AACCTCACTGACATCAAACCTCAACAGTTCAGGAGTTTGCTCAGCTCTCTCTAGACTGGGAGGC
TTCTGGAGGAGGCGTGGCAGCCAGTTTCAGGCTGGCCCGCTGTCCAGCTGGCGGGTTCATGCGGGG
CTGGTCAGCCACAGTGCCGTGAGGCCCAAGGGGCTCTGGTGAGAGGATTATCCACACCCCTC
TACAGTGCCAGAAATCATGACTACGACGTGCCCTCCTGAGGCTCCAGACCGCTCAACTTCTCAGAC
ACTGTGGGCGTGTGTGCTGCGGCCAAGGAACAGCATTTTCCGAAGGGCTCGCGGTGCTGGGTGCT
GGCTGGGGCCACCCACCCTAGCCATACTTACAGCTCGGATATGCTCCAGGACACGGTGGTGGCCTTG
TTCAGCACTCAGCTCTGCAACAGCTCTTGCGGTACAGCGGAGCCCTCACCCCCGATGCTTTGCGCT
GGCTACCTGGACGGAAGGGCTGATGCATGCCAGGGAGATAGCGGGGGCCCTAGTGTGCCAGATGGG
GACACATGGCGCTAGTGGGGTGGTCACTGGGGCGTGGCTGCGCAGAGCCCAATCACCCAGGTGTC
TACGCCAAGGTAGCTGAGTTTCTGGACTGGATCCATGACACTGCTCAGGACTCCCTCCTTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: SgfI-MluI



[View online »](#)

Plasmid Map:


ACCN: NM_001288752

Insert Size: 1167 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001288752.1](#)

RefSeq Size: 2053 bp

RefSeq ORF: 1167 bp

Locus ID: 80975

UniProt ID: [Q9H3S3](#)

Cytogenetics: 11q23.2

Protein Families: Druggable Genome, Protease, Transmembrane

MW: 42.3 kDa

Gene Summary:

This gene encodes a protein that belongs to the serine protease family. Serine proteases are known to be involved in many physiological and pathological processes. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jan 2014]

Transcript Variant: This variant (5) lacks two consecutive in-frame exons in the central coding region, compared to variant 1. The encoded isoform (5) is shorter, compared to isoform 1.