

## Product datasheet for **SC335810**

### VRK2 (NM\_001288839) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	VRK2 (NM_001288839) Human Untagged Clone
Tag:	Tag Free
Symbol:	VRK2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC335810 representing NM_001288839. Blue=Insert sequence Red=Cloning site Green=Tag(s)


```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGTAATGAAAAGACTAGGAATAGATTTACAGAAGATCTCAGGCCAGAATGGTACCTTTAAAAAGTCA
ACTGTCCTGCAATTAGGTATCCGAATGTTGGATGTACTGGAATATATACATGAAAATGAATATGTTTCAT
GGTGATATAAAAGCAGCAAACTACTTTTGGGTTACAAAAATCCAGACCAGGTTTATCTTGACAGATTAT
GGACTTTCCTACAGATATTGTCCAATGGGAACCACAAACAGTATCAGGAAAATCCTAGAAAAGGCCAT
AATGGGACAATAGAGTTTACCAGCTTGGATGCCACAAGGGAGTAGCCTTGTCAGACGAAGTGACGTT
GAGATCCTCGGCTACTGCATGCTGCGGTGGTTGTGTGGGAAACTCCCTGGGAACAGAACCTGAAGGAC
CCTGTGGCTGTGCAGACTGCTAAAACAACTCTGTTGGACGAGCTCCCCCAGTCAGTGCTTAAATGGGCT
CCTTCTGGAAGCAGTTGCTGTGAAATAGCCCAATTTTTGGTATGTGCTCATAGTTTAGCATATGATGAA
AAGCCAACTATCAAGCCCTCAAGAAAATTTGAACCCTCATGGAATACCTTTAGGACCACTGGACTTT
TCCACAAAAGGACAGAGTATAAATGTCCATACTCAAACAGTCAAAAAGTTGATTCACAAAAGGCTGCA
ACAAAGCAAGTCAACAAGGCACACAATAGGTTAATCGAAAAAAAAGTCCACAGTGAGAGAAGCGCTGAG
TCCTGTGCAACATGGAAGTGCAGAAAGAGGAGAACTGATTGGATTGATGAACAATGAAGCAGCTCAG
GAAAGCACAAGGAGAAGACAGAAATATCAAGAGTCTCAAGAACCTTTGAATGAAGTAAACAGTTTCCCA
CAAAAATCAGCTATACACAATCCCAAATCATTATGAGCCTCATCAAGATTTTACCAGTCCAGAT
ATATTCAGAAGTCAAGATCTCCATCTTGGTATAAATACACTCCACAGTCAGCAGGGGATCACAGAC
TTAGAAAGTTCAACTGGACTTTGGCCTACAATTTCCAGTTTACTCTTAGTGAAGAGACAAACGCAGAT
GTTTATTATTATCGCATCATACCTGTCTTTTGTAGTATTCTTGCTTATTTTTTCTCTGA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGCCCGGC
```

Restriction Sites: SgfI-MluI



[View online »](#)

**Plasmid Map:**


**ACCN:** NM\_001288839

**Insert Size:** 1173 bp

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_001288839.1](#)

**RefSeq Size:** 1916 bp

**RefSeq ORF:** 1173 bp

**Locus ID:** 7444

**UniProt ID:** [Q86Y07](#)

**Cytogenetics:** 2p16.1

**Protein Families:** Druggable Genome, Protein Kinase, Transmembrane

**MW:** 44.4 kDa

**Gene Summary:**

This gene encodes a member of the vaccinia-related kinase (VRK) family of serine/threonine protein kinases. The encoded protein acts as an effector of signaling pathways that regulate apoptosis and tumor cell growth. Variants in this gene have been associated with schizophrenia. Alternative splicing results in multiple transcript variants that differ in their subcellular localization and biological activity. [provided by RefSeq, Jan 2014]

**Transcript Variant:** This variant (7) lacks an exon in the 5' region, compared to variant 1, and initiates translation at a downstream in-frame start codon. The encoded isoform (4) has a shorter N-terminus than isoform 1. Variants 6 and 7 encode the same isoform (4).  
**Sequence Note:** This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.