

Product datasheet for SC335813

NPSR1 (NM_001300935) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	NPSR1 (NM_001300935) Human Untagged Clone
Tag:	Tag Free
Symbol:	NPSR1
Synonyms:	ASRT2; GPR154; GPRA; NPSR; PGR14; VRR1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC335813 representing NM_001300935. Blue=Insert sequence Red=Cloning site Green=Tag(s)

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGCCAGCAACTTCACAGAGGGCAGCTTCGATTCCAGTGGGACCGGCAGACGCTGGATTCTTCCCA
GTGGCTTGCACTGAAACAGTGACTTTTACTGAAGTGGTGAAGGAAAGGAATGGGGTTCTTCTACTAC
TCCTTTAAGACTGAGCAATTGATAACTCTGTGGTCTCTTTGTTTTACCATTGTTGGAACTCCGTT
GTGCTTTTTCCACATGGAGGAGAAAGAAGAAGTCAAGAATGACCTTCTTTGTGACTCAGCTGGCCATC
ACAGATCTTTCCACAGGACTGGTCAACATCTTGACAGATATTAATTGGCGATTCACTGGAGACTTCAG
GCACCTGACCTGGTTTGGCGAGTGGTCCGCTATTTGCAGGTTGTGCTGCTCTACGCTCTACCTACGTC
CTGGTGTCCCTCAGCATAGACAGATACCATGCCATCGTCTACCCCATGAAGTTCTTCAAGGAGAAAAG
CAAGCCAGGGTCTCATTGTGATCGCCTGGAGCCTGTCTTTTCTGTTCTCCATTCCACCCTGATCATA
TTTGGGAAGAGGACACTGTCCAACGGTGAAGTGCAGTGTGGCCCTGTGGCCTGACGACTCCTACTGG
ACCCCATACATGACCATCGTGGCCTTCTGGTGTACTTCATCCCTCTGACAATCATCAGCATCATGTAT
GGCATTGTGATCCGAACTATTTGGATAAAAGCAAAACCTACGAAACAGTGATTTCCAACGTCTCAGAT
GGGAACTGTGCAGCAGCTATAACCGAGGACTCATCTCAAAGGCAAAAATCAAGGCTATCAAGTATAGC
ATCATCATATTCTTGCCTTCTCTGCTGTTGGAGTCCATACTTCTGTTGACATTTGGACAATTTT
AACCTCCTTCCAGACACCCAGGAGCGTTTCTATGCCTCTGTGATCATTGACAACCTGCCAGCATTGAAT
AGTGCCATCAACCCCTCATCTACTGTGTCTTCCAGCAGCTCCATCTTTCCCTGCAGGGCTAATAGC
AGTGTCTACCTGTTGGCCTGTGATGTCTCTGCTCTGGGCTCTGGTGTGCTGACCAGTGAGAAGGAGAGC
TGTGAGTCATGGAGAAGGAAGGGAGCAAAGATCACAGGATTCCAGAATGACGTTCCGGGAGAGA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC

Restriction Sites: Sgfl-MluI



[View online >](#)

ACCN:	NM_001300935
Insert Size:	1173 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001300935.1
RefSeq Size:	1713 bp
RefSeq ORF:	1173 bp
Locus ID:	387129
UniProt ID:	Q6W5P4
Cytogenetics:	7p14.3
Protein Families:	Druggable Genome, Transmembrane
MW:	44.3 kDa
Gene Summary:	<p>This gene encodes a member of the vasopressin/oxytocin subfamily of G protein-coupled receptors. The encoded membrane protein acts as a receptor for neuropeptide S and affects a variety of cellular processes through its signaling. Increased expression of this gene in ciliated cells of the respiratory epithelium and in bronchial smooth muscle cells is associated with asthma. Polymorphisms in this gene have also been associated with asthma susceptibility, panic disorders, inflammatory bowel disease, and rheumatoid arthritis. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jul 2014]</p> <p>Transcript Variant: This variant (5) includes an alternate exon in the 3' coding region, which results in a frameshift, compared to variant 1. The encoded isoform (G) has a longer and distinct C-terminus, compared to isoform A.</p>