

Product datasheet for **SC335846**

HARS1 (NM_001289093) Human Untagged Clone

Product data:


Product Type:	Expression Plasmids
Product Name:	HARS1 (NM_001289093) Human Untagged Clone
Tag:	Tag Free
Symbol:	HARS1
Synonyms:	CMT2W; HARS; HRS; USH3B
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC335846 representing NM_001289093. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTGTAGTAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGGCAGAGCGTGCGGCGCTGGAGGAGCTGGTAAACTTCAGGGAGAGCGCGTGCGAGGCCCTCAAGCAG
CAGAAGGCCAGCGCCGAGCTGATCGAGGAGGAGGTGGCGAAACTCCTGAAACTGAAGGCACAGCTGGGT
CCTGATGAAAGCAAACAGAAATTTGTGCTCAAACCCCAAGGATTTTGACATTGCTGGGAACTTTGAT
CCCATGATCCCTGATGCAGAGTGCCTGAAGATCATGTGCGAGATCCTGAGTTCACATTCAGATAGGCGAC
TTCTGGTCAAGGTAACGATCGACGCATTCTAGATGGGATGTTTGTATCTGTGGTGTCTGACAGC
AAGTTCGTACCATCTGCTCCTCAGTAGACAAGCTGGACAAGGTGTCCTGGGAAGAGGTGAAGAATGAG
ATGGTGGGAGAGAAGGGCCTTGACCTGAGGTGGCTGACCGCATTGGGGACTATGTCCAGCAACATGGT
GGGGTATCCCTGGTGAACAGCTGCTCCAGGATCCTAAACTATCCAAAAACAAGCAGGCCTTGAGGGC
CTGGGAGACCTGAAGTTGCTCTTTGAGTACCTGACCTATTTGGCATTGATGACAAAATCTCCTTTGAC
CTGAGCCTTGCTCGAGGGCTGGATTACTACACTGGGGTGTCTATGAGGCAGTGTCTACAGACCCCA
GCCAGGCAGGGGAAGAGCCCTGGGTGTGGGCAGTGTGGCTGCTGGAGGACGCTATGATGGGCTAGTG
GGCATGTTGACCCCAAGGGCGCAAGGTGCCATGTGTGGGCTCAGCATTGGGGTGGAGCGGATTTTC
TCCATCGTGGAACAGAGACTAGAGGCTTTGGAGGAGAAGATACGGACCACGGAGACACAGGTGCTTGTG
GCATCTGCACAGAAGAAGCTGCTAGAGGAAAGACTAAAGCTTGTCTCAGAACTGTGGGATGCTGGGATC
AAGGCTGAGCTGCTGTACAAGAAGAACCCTAAAGCTACTGAACAGTTACAGTACTGTGAGGAGGCAGGC
ATCCCACTGGTGGCTATCATCGGCGAGCAGGAACTCAAGGATGGGGTTCATCAAGTCCGTTCAAGTACG
AGCAGGGAAGAGGTGGATGTCCGAAGAGAAGACCTTGTGGAGGAAATCAAAGGAGAACAGGCCAGCCC
CTCTGCATCTGCTGA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: Sgfl-Mlul



[View online >](#)

Plasmid Map:


ACCN: NM_001289093

Insert Size: 1188 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001289093.1](#)

RefSeq Size: 1964 bp

RefSeq ORF: 1188 bp

Locus ID: 3035

Cytogenetics: 5q31.3

Protein Pathways: Aminoacyl-tRNA biosynthesis

MW: 43.9 kDa

Gene Summary:

Aminoacyl-tRNA synthetases are a class of enzymes that charge tRNAs with their cognate amino acids. The protein encoded by this gene is a cytoplasmic enzyme which belongs to the class II family of aminoacyl-tRNA synthetases. The enzyme is responsible for the synthesis of histidyl-transfer RNA, which is essential for the incorporation of histidine into proteins. The gene is located in a head-to-head orientation with HARSL on chromosome five, where the homologous genes share a bidirectional promoter. The gene product is a frequent target of autoantibodies in the human autoimmune disease polymyositis/dermatomyositis. Several transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Apr 2012]

Transcript Variant: This variant (6) lacks three alternate in-frame exons in the 5' coding region, compared to variant 1. The encoded protein (isoform 6) is shorter, compared to isoform 1.