

Product datasheet for **SC335852**

SGCE (NM_001301139) Human Untagged Clone

Product data:


Product Type:	Expression Plasmids
Product Name:	SGCE (NM_001301139) Human Untagged Clone
Tag:	Tag Free
Symbol:	SGCE
Synonyms:	DYT11; epsilon-SG; ESG
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC335852 representing NM_001301139. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGCAATTGCCCGGTGGTGGGAGCTGGGAGACCCTGTGCTTGGACGGGACAGGGTCGGGGGACACGC
AGGATGAGCCCCGCGACCACTGGCACAATCTTGTCTGACAGGCGAGATTAGTAATGATCCCATAACATTT
AATACAAATTTAATGGGTACCCAGACCGACTGGATGGCTTCGATATATCCAAAGGACACCATATAGT
GATGGAGTCTATATGGGTCCCCAACAGCTGAAAATGTGGGAAGCCAACAATCATTGAGATAACTGCC
TACAACAGGCGCACCTTTGAGACTGCAAGGCATAATTTGATAATTAATATAATGTCTGCAGAAGACTTC
CCGTTGCCATATCAAGCAGAATTCCTTAAGAATATGAATGTAGAAGAAATGTTGGCCAGTGAGGTT
CTTGGAGACTTTCTTGGCGCAGTGAAAAATGTGTGGCAGCCAGAGCGCCTGAACGCCATAAACATCACA
TCGGCCCTAGACAGGGGTGGCAGGGTGCCACTTCCCATAATGACCTGAAGGAGGGCGTTTATGTATG
GTTGGTGCAGATGTCCTGTTTCTTCTTGTTCACGAGAAGTTGAAAATCCACAGAATCAATTGAGATGT
AGTCAAGAAATGGAGCCTGTAATAACATGTGATAAAAAATTTTCGTAATTTTACATTGACTGGTGC
AAAATTTTCATTGGTTGATAAAACAAAGCAAGTGTCCACCTATCAGGAAGTGATTCGTGGAGAGGGGATT
TTACCTGATGGTGGAGAATACAAACCCCTTCTGATTCTTTGAAAAGCAGAGACTATTACACGGATTTT
CTAATTACACTGGCTGTGCCCTCGGCAGTGGCACTGGTCTTTTTCTAATACTTGCTTATATCATGTGC
TGCCGACGGGAAGCGTGGAAAAGAGAAACATGCAAACACCAGACATCCAAGTCCATCACAGTGT
ATTCAGAAATCTACCAAGGAGCTTCGAGACATGTCCAAGAATAGAGAGATAGCATGGCCCTGTCAACG
CTTCTGTGTTCCACCCTGTGACTGGGAAATCATACCTCTTTACACACAGACAATATGATAGCACA
AACATGCCATTGATGCAAACGCAGCAGAACTTGCCACATCAGACTCAGATTCCCCAACAGCAGACTACA
GGTAAATGGTATCCCTGA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: Sgfl-MluI



[View online >](#)

Plasmid Map:


ACCN: NM_001301139

Insert Size: 1191 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001301139.1](#)

RefSeq Size: 1586 bp

RefSeq ORF: 1191 bp

Locus ID: 8910

UniProt ID: [O43556](#)

Cytogenetics: 7q21.3

Protein Families: Transmembrane

MW: 45.1 kDa

Gene Summary:

This gene encodes the epsilon member of the sarcoglycan family. Sarcoglycans are transmembrane proteins that are components of the dystrophin-glycoprotein complex, which link the actin cytoskeleton to the extracellular matrix. Unlike other family members which are predominantly expressed in striated muscle, the epsilon sarcoglycan is more broadly expressed. Mutations in this gene are associated with myoclonus-dystonia syndrome. This gene is imprinted, with preferential expression from the paternal allele. Alternatively spliced transcript variants encoding different isoforms have been found for this gene. A pseudogene associated with this gene is located on chromosome 2. [provided by RefSeq, Oct 2016]

Transcript Variant: This variant (4) lacks two consecutive in-frame exons in the 5' coding region, compared to variant 5. It encodes isoform 4, which lacks an internal in-frame segment and is shorter, compared to isoform 5.