

Product datasheet for **SC335854**

RHD (NM_001282869) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	RHD (NM_001282869) Human Untagged Clone
Tag:	Tag Free
Symbol:	RHD
Synonyms:	CD240D; DIIIc; RH; Rh4; RH30; RHCED; RhDCw; RHDel; RHDVA(TT); RhII; RhK562-II; RhPI; RHPII; RHXIII
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC335854 representing NM_001282869. Blue=Insert sequence Red=Cloning site Green=Tag(s)


```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGAGCTCTAAGTACCCGCGGTCTGTCCGGCGCTGCCTGCCCTCTGGGCCCTAACACTGGAAGCAGCT
CTCATTCTCCTCTTCTATTTTTTTTACCCACTATGACGCTTCCTTAGAGGATCAAAAGGGGCTCGTGGCA
TCCTATCAAGTTGGCCAAGATCTGACCGTGATGGCGGCCATTGGCTTGGGCTTCCTCACCTCGAGTTTC
CGGAGACACAGCTGGAGCAGTGTGGCCTTCAACCTCTTCATGCTGGCGCTTGGTGTGCAGTGGGCAATC
CTGCTGGACGGCTTCTGAGCCAGTTCCTTCTGGGAAGGTGGTCATCACACTGTTCAAGTATTGGCTG
GCCACCATGAGTCTTTGTCGGTCTGATCTCAGTGGATGCTGTCTTGGGAAGGTCAACTTGGCGCAG
TTGGTGGTGGTGGTGGTGGAGGTGACAGCTTTAGGCAACCTGAGGATGGTCATCAGTAATATCTTC
AACACAGACTACCACATGAACATGATGCACATCTACGTGTTCCGAGCCTATTTTGGGCTGTCTGTGGCC
TGGTGCCTGCCAAAGCCTTACCCGAGGGAACGGAGGATAAAGATCAGACAGCAACGATACCCAGTTTG
TCTGCCATGCTGGGCGCCCTCTTCTGTGGATGTTCTGGCCAAGTTTCAACTCTGCTCTGCTGAGAAGT
CCAATCGAAAGGAAGAATGCCGTGTTCAACACCTACTATGCTGTAGCAGTCAGCGTGGTGACAGCCATC
TCAGGGTCATCCTTGGCTCACCCCAAGGGAAGATCAGCAAGACTTATGTGCACAGTGGGTGTTGGCA
GGAGGCGTGGCTGTGGGTACCTCGTGTACCTGATCCCTTCTCCGTGGCTTGGCATGGTGTGGGCTT
GTGGCTGGGCTGATCTCCGTCCGGGGAGCCAAGTACCTGCCGGGTGTTGTAACCGAGTGTGGGGATT
CCCCACAGCTCCATCATGGGCTACAACCTCAGCTTGTGGGTCTGCTTGGAGAGATCATCTACATTGTG
CTGCTGGTGTGATACCGTCGGAGCCGCAATGGCATTTCCTCATTTGGCTGTTGGATTTAAGCAA
AAGCATCAAGAAAAACAAGGCCTGTTCAAAAACAAGCAACTTCTCTACTGTTGCCTGCATTGTGA
CGTGAGAAACGCTCATGA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```



[View online >](#)

Restriction Sites: Sgfl-MluI

Plasmid Map:



ACCN: NM_001282869

Insert Size: 1191 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001282869.1](#)

RefSeq Size: 2683 bp

RefSeq ORF: 1191 bp

Locus ID: 6007

UniProt ID: [Q02161](#)

Cytogenetics: 1p36.11

Protein Families: Transmembrane

MW: 43.3 kDa

Gene Summary:

The Rh blood group system is the second most clinically significant of the blood groups, second only to ABO. It is also the most polymorphic of the blood groups, with variations due to deletions, gene conversions, and missense mutations. The Rh blood group includes this gene, which encodes the RhD protein, and a second gene that encodes both the RhC and RhE antigens on a single polypeptide. The two genes, and a third unrelated gene, are found in a cluster on chromosome 1. The classification of Rh-positive and Rh-negative individuals is determined by the presence or absence of the highly immunogenic RhD protein on the surface of erythrocytes. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2008]

Transcript Variant: This variant (5) lacks two alternate exons compared to variant 1, that causes a frameshift. The resulting isoform (4) has a shorter and distinct C-terminus compared to isoform 1. Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.