

Product datasheet for **SC335859**

Bestrophin 3 (BEST3) (NM_001282614) Human Untagged Clone

Product data:


Product Type:	Expression Plasmids
Product Name:	Bestrophin 3 (BEST3) (NM_001282614) Human Untagged Clone
Tag:	Tag Free
Symbol:	BEST3
Synonyms:	VMD2L3
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC335859 representing NM_001282614. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTGTAGTGAACCGTCAGAATTTTGTAAATACGACTCACTATAGGGCGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGACTGTCACCTACTCCAGTAAAGTAGCAAAATGCAACTTTTTTTGGATTCATAGGTTACTCCTCAAG
TGGAGAGGCAGCATCTACAACTACTGTACAGGAATTTATTGTTTTGCTGTTCTTTATACAGCAATA
AGTTTGGTATACAGATTGTTACTTACAGGAGTCCAAAACGTTACTTTGAAAAATATCAATTTACTGT
GACAGATATGCTGAACAAATCCAGTAACCTTTGTGCTGGGTTTTATGTTACTCTGGTAGTGAACCGA
TGGTGAACCCAGTTTGTGAATTTGCCCTGGCCAGACAGGCTAATGTTCCATCTCTAGCAGTGTTCAC
GGAAGCGACGAGCACGGGCGCTGCTTAGAAGGACGCTGATGCGCTACGTCATCTCACCTCCCTGCTC
ATCTTTCGCTCGGTGAGCACTGCTGTGTACAAAAGATTTCCACAATGGACCACGTGTTGAAGCAGGT
TTTATGACAACAGATGAAAGGAAATTTCAACCACCTCAAGTCTCCTCATCTGAAATATTGGGTTCCA
TTCATCTGGTTTTGAAATCTTGCAACTAAAGCCCGGAATGAAGGTAGAATCAGAGACAGTGTGATCTG
CAATCATTGATGACTGAAATGAATCGATACCGCTCTTGGTGCAGCCTCTTATTGCGTTATGACTGGGT
GGGATCCGCTGGTTTACACCCAGGTTGCTACTTGTGCTGTATACCTCTTCTTTGCGTGCCTGATT
GGACGCCAGTTTTTGGATCCCACCAAGGCTACGCAGGGCATGACTTGGATCTTTACATTTCCCATCTC
ACCCTCTACAATTCTTCTTCTATGCAGGATGGCTTAAGGTAGCAGAGCAGCTTATCAACCTTTTGGG
GAAGATGATGATGATTTTGAACCTAACTGGTGCATTGACAGAAATTTGCAGGCTCTCTTTTAGCTGTG
GACGAAATGCACATGAGCTTACCCAAGATGAAGAAGGACATTTACTGGGACGATTCTGCTCGCCCA
CCATACACATTGGCAGCTGCTGACTACTGCATACCCTCATTTCTGGGGTCAACAGTCCAGATGGGAAA
CAGATGCCTAAGAATGAGTGGAAGATGGAAGATATAAAAAATCCCACTCCCTCAGCCTCAATTCCAGTGT
GCTAAGAGTGATCCTGGTGGTTAA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
```

Restriction Sites: Sgfl-MluI



[View online »](#)

Plasmid Map:


ACCN: NM_001282614

Insert Size: 1197 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001282614.1](#)

RefSeq Size: 1607 bp

RefSeq ORF: 1197 bp

Locus ID: 144453

UniProt ID: [Q8N1M1](#)

Cytogenetics: 12q15

Protein Families: Ion Channels: Other, Transmembrane

MW: 46.6 kDa

Gene Summary:

BEST3 belongs to the bestrophin family of anion channels, which includes BEST1 (MIM 607854), the gene mutant in vitelliform macular dystrophy (VMD; MIM 153700), and 2 other BEST1-like genes, BEST2 (MIM 607335) and BEST4 (MIM 607336). Bestrophins are transmembrane (TM) proteins that share a homology region containing a high content of aromatic residues, including an invariant arg-phe-pro (RFP) motif. The bestrophin genes share a conserved gene structure, with almost identical sizes of the 8 RFP-TM domain-encoding exons and highly conserved exon-intron boundaries. Each of the 4 bestrophin genes has a unique 3-prime end of variable length (Stohr et al., 2002 [PubMed 12032738]; Tsunenari et al., 2003 [PubMed 12907679]).[supplied by OMIM, Mar 2008]

Transcript Variant: This variant (4) contains an alternate 3' terminal exon and it thus differs in the 3' coding region and 3' UTR, compared to variant 1. The encoded isoform (4) has a distinct and shorter C-terminus, compared to isoform 1.