

## Product datasheet for **SC335862**

### ZNF302 (NM\_001289186) Human Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** ZNF302 (NM\_001289186) Human Untagged Clone  
**Tag:** Tag Free  
**Symbol:** ZNF302  
**Synonyms:** HSD16; MST154; MSTP154; ZNF135L; ZNF140L; ZNF327  
**Vector:** pCMV6-Entry (PS100001)  
**Fully Sequenced ORF:** >SC335862 representing NM\_001289186.  
Blue=Insert sequence Red=Cloning site Green=Tag(s)


```
ATGTCAGGTCACATTTAGTGATGTGGCTATAGACTTCTCTCATGAAGAGTGGGCATGCCTAGATTCT
GCTCAGAGGGACTTATACAAGGATGTGATGGTCCAGAATTATGAGAACCTGGTCTCTGTAGGTCTTTCC
GTAACAAGCCATATGTGATCATGTTATTGGAGGATGGAAAAGAGCCCTGGATGATGGAGAAAAAAGT
TCAAAAAGATTGGGAATCAAGATGGGAAAACAAGGAATTATCAACAAGAAGGATATTTATGATGAAGAT
TCACCCCAACCAGTAACAATGGAAAAAGTTGTAACAACAAGTTATGAATTTTCAAATTCTAATAAGAA
TTGGAATATACAGAATGCGACACATTTAGAAGCACCTTTTATTCAAAGTCTACTCTTTCTGAACCACAA
AACAAATTCGCTGAAGGGAATTCACACAAATATGATATATTAAAGAAGAATTTATCAAAAAAGTCAGTT
ATAAAAAGTGAGAGAATAAATGGTGGAAAAGAACTTTTAAATTCTAATAAAAAGTGGGGCAGCCTTCAAC
CAGAGCAAATCTCTTACCCTTCCCAGACTTGTAAATAGAGAGAAAATCTATACATGCAGTGAATGTGGG
AAAGCCTTTGGCAAACAGTCAATCCTCAGTCGCCACTGGAGAATTCATACAGGAGAGAAGCCCTATGAA
TGTCGTGAATGTGGGAAGACTTTTAGCCATGGTTCATCCCTTACACGACATCAGATAAGCCATAGTGA
GAGAAACCTTACAAATGCATTGAATGTGGGAAGGCCTTTAGCCATGGCTCATCACTTACTAACCATCAG
AGCACTCACACGGGAGAGAAAACCGTATGAATGTATGAACTGTGGAAGTCTTTTAGTTCGTGTGTCCCTT
CTCATTACAGCATCTAAGAATTCATACGCAAGAAAACGCTATGAGTGTGATATGTGGAAGGCCTTC
ATTCATAGTTTCGTCTCTCATTACCATCAGAAAAGCCATACTGGAGAGAAGCCTTATGAATGTAGAGAA
TGTGGGAAAGCTTTCTGCTGTAGCTCACACCTTACTCAACATCAAAGAATTCACAGTATGAAGAAAAA
TATGAATGCAACAAATGTCTCAAGGTCTTTAGTAGCTTCTCATTCTTGTCAACATCAGAGTATTCAT
ACTGAAGAAAACCGTTTGAAGTTAG
```

**Restriction Sites:** SgfI-MluI



[View online »](#)

## Plasmid Map:



ACCN: NM\_001289186

Insert Size: 1200 bp

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM\\_001289186.1](#)

RefSeq Size: 2586 bp

RefSeq ORF: 1200 bp

**Locus ID:** 55900

**UniProt ID:** [Q9NR11](#)

**Cytogenetics:** 19q13.11

**Protein Families:** Transcription Factors

**MW:** 45.9 kDa

**Gene Summary:** This gene encodes a member of the zinc-finger protein family. The encoded protein contains seven C2H2-type zinc fingers and a KRAB domain, but its function has yet to be determined. Alternatively spliced transcript variants have been described. [provided by RefSeq, Mar 2014]  
Transcript Variant: This variant (8) differs in the 5' UTR and lacks an alternate in-frame exon compared to variant 3. The resulting isoform (c) has the same N- and C-termini but is shorter compared to isoform a. Variants 1, 2, 5, 8, and 9 all encode the same isoform (c).