

Product datasheet for SC335879

DOC2A (NM_001282063) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	DOC2A (NM_001282063) Human Untagged Clone
Tag:	Tag Free
Symbol:	DOC2A
Synonyms:	Doc2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC335879 representing NM_001282063. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```


GTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGAATTCGTGCGACTGGATCCGGTA
CCGAGGAGATCTGCCGCCGCGATCGCCGGCGCGCC
ATGAGGGGCCGAGGGGCGATCGCATGACCATCAACATCCAGGAGCACATGGCCATCAACGTGTGCCCC
GGGCCCATCCGGCCCATCCGCCAGATCTCTGACTACTTCCCCGGGGACCAGGACCTGAAGGGGGCGGC
GGGGGGCGGGGAGGCCCGCCCATCTGGTCCCCCTGGCTCTGGCCCCCTGCAGCCCTCCTTGGG
GCCACCACGCCTGAGGATGGTGGGAGGTGGACAGCTATGACTCGGATGATGCCACCGCCCTAGGCAGC
CTGGAGTTTGACCTTCTCTACGACCGGGCTCCTGCACCTGCACTGTAGCATCCTCAGGCCAAGGGC
CTCAAGCCATGGATTTCAATGGCTCGCCGACCCTACGTCAAGCTGCACCTGCTGCCTGGAGCCTGT
AAGGCCAATAAGCTAAAAACGAAGACTCAGAGGAACACTGAATCCCGTGTGGAATGAGGACCTGACT
TACAGCGGGATCACAGATGACGACATCACGCACAAGGTGCTCAGGATCGCCGTCTGTGATGAGGACAAG
CTGAGTCAACAATGAGTTTATTGGGGAGATCCGCGTGCCCTCCGCGCCCTCAAGCCTTCGCAGAAGAAG
CATTTTAAACATCTGCCTCGAGCGCCAGGTCCCCTGGCGTCCCCCTTTCCATGTGAGCGGCGCTGAGG
GGCATCTCCTGTTATCTGAAGGAGTTGGAGCAGGCGGAGCAGGGGAGGGGCTGCTGGAGGAGCGTGGC
CGCATCCTGCTGAGTCTCAGCTACAGCTCGCGGCCCGGGGACTGCTGGTAGGCATCTTGCCTGCGCC
CATCTGGCTGCCATGGACGTCAACGTTACTCGGACCCCTACGTCAAGACGTACCTGAGGCCGATGTG
GACAAGAAATCCAAGCATAAGACGTGTGTGAAGAAGAAGACTCAACCCAGAATTTAACGAGGAGTTT
TTCTACGAGATAGAGCTCTCCACTCTGGCCACCAAGACCCTGGAAGTCAACGCTCTGGGACTATGACATT
GGCAAATCCAATGACTTCATTGGTGGCGTGTCCCTGGGGCCAGGTGCCCGAGGCGAGGCTCGGAAGCAC
TGGAGTACTGCCTGCAGCAGCCGACGACCCCTGGAGCGCTGGCACACCCTGACCAGTGAAGTGCCT
CCTGCGGCCGGGGCTCTGTCTCAGCCTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC

```

Restriction Sites: Ascl-MluI



[View online »](#)

Plasmid Map:


ACCN: NM_001282063

Insert Size: 1203 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001282063.1](#)

RefSeq Size: 2090 bp

RefSeq ORF: 1203 bp

Locus ID: 8448

UniProt ID: [Q14183](#)

Cytogenetics: 16p11.2

Protein Families: Secreted Protein

MW: 44 kDa

Gene Summary:

There are at least two protein isoforms of the Double C2 protein, namely alpha (DOC2A) and beta (DOC2B), which contain two C2-like domains. DOC2A and DOC2B are encoded by different genes; these genes are at times confused with the unrelated DAB2 gene which was initially named DOC-2. DOC2A is mainly expressed in brain and is suggested to be involved in Ca(2+)-dependent neurotransmitter release. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Aug 2013]

Transcript Variant: This variant (2) contains an alternate terminal exon in the 5' UTR compared to variant 1. Variants 1 through 4 encode the same protein.