

Product datasheet for **SC335895**

ATE1 (NM_001288734) Human Untagged Clone

Product data:


Product Type:	Expression Plasmids
Product Name:	ATE1 (NM_001288734) Human Untagged Clone
Tag:	Tag Free
Symbol:	ATE1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC335895 representing NM_001288734. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGATTCCACAATGGATGATGCTGTTGCGGGTGACTTTGCATTGATAAAATAAACTGGATATACAGTGT
GATCTTAAACACTCAGTGATGACATCAAAGAGAGTTTAGAGAGTGAAGGAAAAAATTCAAAGAAAGAA
GAACCTCAGGAATTACTTCAGTCACAAGATTCGTAGGAGAGAAGTTGGGCTCTGGTGAACCGTCACAT
TCAGTTAAAGTTCACACAGTTCCTAAGCCAGGCAAAGGGGCTGATTTGAGTAAGCCTCCATGTCGAAAA
GCAAAGGAAATCCGAAAGAAAGGAAAGGTTAAACTAATGCAGCAGAACCCAGCTGGAGAACTTGAG
GGTTTCCAGGCTCAAGGTCACCCACCATCTTTGTTTCCACCAAAGGCTAAATCCAACCGCCAAAATCA
CTCGAAGATTTAATTTTGGAGTCTTTACCAGAGAATGCATCACACAAGTTAGAGGTGAGGGTGGTGAGA
TCATCTCCACCAAGTTCGCAGTTCAAAGCCACACTTCTGGAGTCTTACCAGGTCTATAAACGTTACCAG
ATGGTTATTCACAAGAACCACCTGATACGCCAACCGAAAGCCAGTTCACAAGATTCCTTTGCAGTTCA
CCCTTGGAGGCAGAGACTCCCCCTAATGGGCCAGATTGTGGCTATGGCTCCTTTCACCAGCAGTACTGG
CTTGACGGAAAGATCATTGCTGTGGGGTGATTGACATCTCCCAAAGTGTATCATCTGTGATTTG
TACTACGATCCTGATTATTCGTTTTGTCTTTGGGCGTCTACTCTGCACTACGAGAAATTGCTTTACT
AGGCAGCTTCATGAGAAAATTCTCAACTCAGCTATTATTATATGGGTTTCTACATTCATCATGTCCC
AAGATGAAATATAAGGGTCAGTATAGACCTTCTGATTTGCTGTGCCCTGAGACATATGTTTGGGTACCC
ATTGAGCAATGCCTGCCTTCACTTGAAAACCTCAAGTACTGCCGTTTCAACCAGGACCCAGAAGCAGTG
GATGAGGATCGCAGTACGGAACCTGACCGATTGCAGGTGTTTCAAGAGAGCCATCATGCCTTACGGT
GTTTATAAGAAACAGCAGAAAGACCCAAGTGAGGAGGCTGCTGTTCTGCAGTACGCCAGCCTGGTGGGG
CAGAAGTGCTCCGAGCGGATGCTGCTGTTCAGAAACTGA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: SgfI-MluI



[View online »](#)

Plasmid Map:


ACCN: NM_001288734

Insert Size: 1212 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001288734.1](#)

RefSeq Size: 4826 bp

RefSeq ORF: 1212 bp

Locus ID: 11101

UniProt ID: [O95260](#)

Cytogenetics: 10q26.13

MW: 45.8 kDa

Gene Summary:

This gene encodes an arginyltransferase, an enzyme that is involved in posttranslational conjugation of arginine to N-terminal aspartate or glutamate residues. Conjugation of arginine to the N-terminal aspartate or glutamate targets proteins for ubiquitin-dependent degradation. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Dec 2013]

Transcript Variant: This variant (3) differs in its 5' UTR and initiates translation at a downstream start codon, compared to variant 1. The encoded isoform (3) has a shorter N-terminus, compared to isoform 1. Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.