

## Product datasheet for SC335898

### MAEL (NM\_001286377) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	MAEL (NM_001286377) Human Untagged Clone
Tag:	Tag Free
Symbol:	MAEL
Synonyms:	CT128; SPATA35
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC335898 representing NM_001286377. Blue=Insert sequence Red=Cloning site Green=Tag(s)


```
GCTCGTTTGTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGCCGAACCGTAAGGCCAGCCGGAATGCTTACTATTTCTTCGTGCAGGAGAAGATCCCCGAACACGG
CGACGAGGCCCTGCCTGTGGCTCGCGTTGCTGATGCCATCCCTTACTGCTCCTCAGACTGGCGAAACCT
GTTTTACACCACTGAGGAGGCCAGGCATGCTTGTACCAAAGCAGAATGTTTCACCTCCAGATATGTCA
GCTTTGTCTTTAAAAGGTGATCAAGCTCTCCTTGGAGGCATTTTTATTTTTTGAACATTTTTAGCCAT
GGCGAGCTACCTCCTCATTGTGAACAGCGCTTCCCTCCCTTGTGAAATTGGCTGTGTTAAGTATTCTC
CAAGAAGTATTATGGCAGATTTCCACAGTTTTATAAATCCTGGTAAATTCACGAGGATTTTCGATTT
CATTGTCAAGGCTGCAAGTGATTCTAGTCACAAGATTCCTATTTCAAATTTGAACGTGGGCATAACCAA
GCAACTGTGTTACAAAACCTTTATAGATTTATTCATCCCAACCCAGGGAAGTGGCCACCTATCTACTGC
AAGTCTGATGATAGAACCAGAGTCAACTGGTGTGTTGAAGCATATGGCAAAGGCATCAGAAATCAGGCAA
GATCTACAACCTTCTCACTGTAGAGGACCTGTAGTGGGGATCTACCAACAAAAATTTCTCAAGGAGCCC
TCTAAGACTTGGATTGAAAGCCTCTAGATGTGGCCATGTGGGATTATTCTAGCAACACAAGGTGCAAG
TGGCATGAAGAAAATGATATTCTTCTGTGCTTAGCTGTTGCAAGAAGATTGCGTACTGCATCAGT
AATTCTCTGGCCACTCTCTTTGGAATCCAGCTCACAGAGGCTCATGTACCACTACAAGATTATGAGGCC
AGCAATAGTGTGACACCCAAAATGGTTGATTGGATGCAGGGCGTTACCAGAAGTAAGGGTTGGGAGT
TCAGGATTCTCTCATTTCAACTTCTAATGAGGAACAAAGATCAAACACACCCATTGGTGACTACCCA
TCTAGGGCAAAAATTTCTGGCCAAAACAGCAGCGTTTCGGGGAAGAGGAATTACCCGCTTACTAGAGAGC
ATTTCCAATTCTCCAGCAATATCCACAAATCTCCAAGTGTGACACTTCACTCTCACCTTACATGTCC
CAAAAAGATGGATACAAATCTTCTCTTCTTATCTTAA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: Sgfl-MluI



[View online »](#)

## Plasmid Map:



ACCN: NM\_001286377

Insert Size: 1212 bp

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM\\_001286377.1](#)

RefSeq Size: 1849 bp

RefSeq ORF: 1212 bp

Locus ID: 84944

UniProt ID: [Q96JY0](#)

Cytogenetics: 1q24.1

MW: 45.6 kDa

**Gene Summary:**

Plays a central role during spermatogenesis by repressing transposable elements and preventing their mobilization, which is essential for the germline integrity. Acts via the piRNA metabolic process, which mediates the repression of transposable elements during meiosis by forming complexes composed of piRNAs and Piwi proteins and governs the methylation and subsequent repression of transposons. Its association with piP-bodies suggests a participation in the secondary piRNAs metabolic process. Required for the localization of germ-cell factors to the meiotic nuage (By similarity).[UniProtKB/Swiss-Prot Function]

Transcript Variant: This variant (2) lacks an in-frame exon in the 5' coding region compared to variant 1. The encoded isoform (2) is shorter than isoform 1. Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.